

Medical Device Details

Sponsor:	Mango Medtech Limited 12 Madden Street Auckland New Zealand
WAND Reference:	211119-WAND-6XZ4OP
Sponsors Own Reference:	AHP300 Transport Ventilator
GMDN:	Transport electric ventilator [36289]
Class:	IIa
Intended Purpose:	The AHP300 is intended to be used as an electrically controlled emergency ventilator that can use an external compressed gas source or its internal compressor. This ventilator is designed to provide emergency respiratory support by means of a face mask or tube inserted into the patient's airway. The ventilator is intended for use on patients weighing greater than 5kg. The ventilator is intended to be used in the environments associated with emergency medical services (EMS), inter-hospital transport and hospital facility usage by qualified, trained personnel under the direction of a physician.
Manufacturer:	Allied Healthcare Products Inc 1720 Sublette Avenue St Louis Missouri United States
ARTG ID:	<ul style="list-style-type: none"> › The device is supplied unsterilised. › The device is intended to be invasive. › The device is not intended for single use. › The device is an active device. › The device does not contain material or ingredients of microbial origin. › The device does not contain material or ingredients of recombinant origin. › The device does not contain material or ingredients manufactured or formulated using a genetically modified organism. › The device does not contain material or ingredients of human origin. › The device does not contain human blood or its components. <p>☐ › The device consists of: System</p> <p>☐ Components:</p> <ul style="list-style-type: none"> › Component: Battery 12v (Reusable), Reusable 12v battery for AHP300 Ventilator (3-year life). › Component: Cuffed Oxygen Mask (Disposable), Disposable cuffed oxygen mask in adult and child variants for AHP300 Ventilator. › Component: Exhalation Filter (Disposable), Disposable exhalation filter for AHP300 Ventilator. › Component: External Battery Pack (Reusable), External battery pack for AHP300 Ventilator. › Component: Replacement Air Hose (Reusable), Reusable air hose (1.8m) with DISS for AHP300 Ventilator. › Component: Replacement Oxygen Hose (Reusable), Reusable replacement Oxygen Hose (1.8m) with DISS. › Component: Single Limb Vent Circuit (Disposable), Disposable Adult or Paediatric Single Limb Vent Circuit (1.8m) for AHP300 Ventilator. › The device does not contain material or ingredients of animal origin. › The device is not medicated. › The device is not formulated. › The product does not contain a medicine that has consent for marketing in New Zealand. › The product does not contain a medical device which incorporates a medicine as an integral part and that has an action ancillary to the device.
Status:	Active