



endoscopic  
visualization,  
reimagined



endoluxe  
WIRELESS VISUALIZATION SYSTEM

# Introducing endoluxe

WIRELESS VISUALIZATION SYSTEM

A **wireless visualization system** with unparalleled functionality and freedom for any analog endoscope.



- The Orb is a 5 mega pixel camera with a touch screen display and proprietary 9000 Lux LED light source



- The Charger conveniently stores and recharges the orb with a reliable docking interface



- The Orb quickly couples to virtually any analog endoscope



- Wirelessly pairs with the Receiver to capture images and video while streaming to your remote display or our cloud application

## Viewer

Orb displays video and transmits wirelessly

Light source magnetically snaps onto scope

Orb couples with rigid or flexible analog endoscopes



## Charger

Accepts images and provides connection to display



## Remote Display



Images and videos upload wirelessly to cloud or download to computer via USB type C port

Endoluxe is the ultimate endoscopy wireless visualization platform with uncompromised quality and cost efficiency. The intuitive operation and simple setup make it the ideal solution for all your visualization needs.

**See the light. Go to [endoluxe.com](http://endoluxe.com)**

