

Implant Company	Model	Platform	Diameter (mm)	SKU No	SmartPeg
NOBEL BIOCARE	Brånemark system	NP		100355	Type 02
	Brånemark system	RP		100353	Type 01
	Brånemark system	WP		100357	Type 03
	Brånemark system Novum			100353	Type 01
	Brånemark system Shorty	RP		100353	Type 01
	Brånemark system Shorty	WP		100357	Type 03
	NobelActive		3	100545	Type 62
	NobelActive	NP	3.5	100543	Type 60
	NobelActive	RP	4.3	100544	Type 61
	NobelActive	RP	5	100544	Type 61
	NobelActive	WP	5.5	100350	Type 04
	NobelDirect	NP		100366	Type 08
	NobelDirect	RP		100368	Type 09
	NobelDirect	WP	3.75	100368	Type 09
	NobelParallel Conical Connection	NP	4.3	100543	Type 60
	NobelParallel Conical Connection	RP	5	100544	Type 61
	NobelParallel Conical Connection	RP	5.5	100544	Type 61
	NobelParallel Conical Connection	WP	3.5	100350	Type 04
	NobelReplace Conical Connection	NP	4.3	100543	Type 60
	NobelReplace Conical Connection	RP	5	100544	Type 61
	NobelReplace Conical Connection	RP		100544	Type 61
	NobelReplace Conical Connection PMC	NP		100543	Type 60
	NobelReplace Conical Connection PMC	RP		100544	Type 61
	NobelReplace Platform Shift	NP		100359	Type 12
	NobelReplace Platform Shift	RP		100361	Type 13
	NobelReplace Platform Shift	WP		100361	Type 13
	NobelReplace Straight Groovy	NP		100359	Type 12
	NobelReplace Straight Groovy	RP		100361	Type 13
	NobelReplace Straight Groovy	WP		100361	Type 13
	NobelReplace Tapered Groovy		6	100462	Type 40
	NobelReplace Tapered Groovy	NP		100359	Type 12
	NobelReplace Tapered Groovy	RP		100361	Type 13
	NobelReplace Tapered Groovy	WP		100361	Type 13
	NobelSpeedy Groovy	NP		100355	Type 02
	NobelSpeedy Groovy	RP		100353	Type 01
	NobelSpeedy Groovy	WP		100357	Type 03
	NobelSpeedy Replace		6	100462	Type 40
	NobelSpeedy Replace	NP		100359	Type 12
	NobelSpeedy Replace	RP		100361	Type 13
	NobelSpeedy Replace	WP		100361	Type 13
	NobelSpeedy Shorty	RP		100353	Type 01

	NobelSpeedy Shorty	WP		100357	Type 03
	Replace Select Straight	NP		100359	Type 12
	Replace Select Straight	RP		100361	Type 13
	Replace Select Straight	WP		100361	Type 13
	Replace Select Tapered		6	100462	Type 40
	Replace Select Tapered	NP		100359	Type 12
	Replace Select Tapered	RP		100361	Type 13
	Replace Select Tapered	WP		100361	Type 13
	Replace Select Tapered PMC		6	100462	Type 40
	Replace Select Tapered PMC	NP		100359	Type 12
	Replace Select Tapered PMC	RP		100361	Type 13
	Replace Select Tapered PMC	WP		100361	Type 13
	Replace Select TC	NP		100359	Type 12
	Replace Select TC	RP		100361	Type 13