




DODO AIR Face Mask

ppe GERMANY

20 Pieces

DODO AIR Gesichtsmasken werden in Berlin (Deutschland) von PPE Germany GmbH produziert. Sie sind zum einmaligen Gebrauch bis zu 8 Stunden bestimmt. Bitte waschen Sie sich vor dem Anlegen der Maske gründlich die Hände. Bitte tauschen Sie die Maske nach Gebrauch sofort gegen eine neue, trockene, saubere und gut sitzende Maske aus. Die Maske ist eine in der Europäischen Union zugelassene FFP2 Maske und entspricht der europäischen DIN EN 149:2001+A1:2009 und der Europäischen PSA-Verordnung EU 2016/425. Haltbarkeit der Maske: 2 Jahre. Mehr Informationen: <https://ppegermany.de/dodo-usa>

DODO AIR face masks are produced in Berlin (Germany) by PPE Germany GmbH. They are intended for single use up to 8 hours. Please wash your hands thoroughly before putting on the mask. Please replace the mask immediately after use with a new, dry, clean and well-fitting mask. The mask is an FFP2 mask approved in the European Union and complies with the European DIN EN 149:2001+A1:2009 and the European PPE Regulation EU 2016/425. Shelf life of the product: 2 years. More information: <https://ppegermany.de/dodo-usa>



MADE IN GERMANY



<https://ppegermany.de/dodo-usa>



Improved Breathability



Adjustable shape



Less fogged glasses



FFP2



DODO AIR

Ready to go

Each mask in this box is individually sealed, making it perfect to take one at a time. Ideal for dispensers, offices, schools and reception areas.

Specifications

Product name: DODO AIR face mask
 EAN: 4260763950424
 Product quantity: 20 pieces per box
 Uses: Dust proof, professional and protective.
 Storage conditions: Temperature from 5 °C to +25 °C and at a maximum humidity of 60%



Manufacturer information

Execution standard: FFP2
 Manufacturer: PPE Germany GmbH, Germany
 Mertsstraße 63-115 · D-13587 Berlin
 Made in Germany – synonymous of product quality, durability and reliability

ppe GERMANY

20 Pieces

DODO AIR Face Mask



DODO AIR Face Mask

ppe GERMANY



DODO AIR

Face Mask

Breathable · Comfortable · Disposable
German engineering and quality assurance tested



20 Pieces



20 Pieces


DODO AIR Face Mask




DODO AIR Face Mask

ppe GERMANY


20 Pieces




1 Place the mask in your hands with the top side facing upwards.



2 Put on the protective mask and press the metal part on the bridge to the shape of your nose.



3 Pull the bottom of the mask slightly down so that the nose, mouth and chin are well covered.



4 Check the fit of the mask so that no secondary air is drawn in. Check facial hair for air leaks. If necessary, correct the fit until the mask fits comfortably and seals correctly.

FFP2

MADE IN GERMANY

For product quality, durability and reliability

21023-16-FFP2-DODO-AIR-medialblack-USA-20ir

