

## Aggrecan Antibody to C-terminal Neoepitope NITEGE Mouse Monoclonal Antibody

**Catalog Number 1042003**

*For research use only*

### INTRODUCTION

Aggrecan is a member of a family of large, aggregating proteoglycans (also including versican, brevican and neurocan) which is found in articular cartilage. Aggrecan is composed of three major domains: G1, G2, and G3. Between the G1 and G2 domains there is an interglobulin region (IGD). The IGD region is the major site of cleavage by specific proteases like metalloproteinases (MMPs) and aggrecanase. Aggrecan cleavage has been associated with a number of degenerative diseases including rheumatoid arthritis and osteoarthritis. There is evidence that this family of proteoglycans modulates cell adhesion, migration, and axonal outgrowth in the CNS.

This aggrecan degradation product usually remains within the tissue still complexed to hyaluronan and link protein. Its release from the cartilage usually signals that there has been extensive catabolism of aggrecan, which allows large complexes containing this metabolite to be released from the tissue.

### IMMUNOGEN

NITEGE synthetic peptide conjugate.

### CLONE

BC-13

### HOST

Mouse

### MYELOMA

x63-Ag8.653

### ISOTYPE

IgG1

### LIGHT CHAIN TYPE

kappa

### SPECIFICITY

Recognizes the aggrecanase (ADAMTS-1, -4 & -5)-generated C-terminal neoepitope ...TEGE after cleavage between amino acids ..EGE373 and 374ARG.. within the interglobular domain of aggrecanase-catabolised aggrecan (Human aggrecan sequence enumeration). This antibody cross-reacts with Human, Cow, Dog, Horse, Pig.

#### International

MD Bioproducts GmbH  
Tel: +41-44 986 2628  
info@mdbioproducts.com



[www.mdbioproducts.com](http://www.mdbioproducts.com)

#### North America

MD Biosciences Bioproducts, LLC  
Tel: 651-789-6535  
products@mdbiosciences.com

<b>PURITY</b>	Affinity purified on protein G
<b>STORAGE BUFFER</b>	PBS, no preservatives
<b>FORM</b>	Liquid
<b>CONCENTRATION</b>	Please see vial label for concentration.
<b>APPLICATION</b>	<ul style="list-style-type: none"><li>• Western-Blotting Suggested dilution: 1:100 Detects a a band of approximately 60-70 kDa</li><li>• ELISA</li><li>• IHC</li></ul>
<b>TECHNICAL NOTES</b>	Samples are usually deglycosylated using 0.01 Units chondroitinase ABC (Sigma), 0.01 Units Keratanase (Seikagaku) and 0.0001 Units Keratanase II (Seikagaku) per 10µg S-GAG of non-deglycosylated aggrecan for optimal epitope recognition in SDS-PAGE and immunohistochemistry (1, 2).
<b>STORAGE</b>	The antibody is stable until the expiry date given on the label if stored at -20 °C. Repeated freezing and thawing should be avoided. Aliquoting is recommended.
<b>EXPIRATION</b>	See Vial Label
<b>REFERENCES</b>	<ol style="list-style-type: none"><li>1. Little CB et al. Matrix metalloproteinases are involved in C-terminal and interglobular domain processing of cartilage aggrecan in late-stage cartilage degradation. Matrix Biol 21:271-88 (2002).</li><li>2. Caterson B et al. Mechanisms involved in cartilage proteoglycan catabolism. Matrix Biol 19:333-44 (2000)</li></ol>

**International**  
MD Bioproducts GmbH  
Tel: +41- 44 986 2628  
info@mdbioproducs.com

   
[www.mdbioproducs.com](http://www.mdbioproducs.com)

**North America**  
MD Biosciences Bioproducts, LLC  
Tel: 651-789-6535  
products@mdbiosciences.com