

Flexible Disc

For Invisible Partial Denture


NOBILCAM Flexible Discs (Nylon resin) are intended for aesthetically superior semi-flexible removable partial denture frameworks, which blend perfectly with the patient's gums, providing better function, aesthetics, comfort and durability.



INVISIBLE

FLEXIBLE

COMFORTABLE



Advantages:

- Esthetic, comfortable and lightweight
- 100% digital for accurate fit
- Occlusal rests reduce tissue-borne pressure to reduce sore spots and offer excellent retention




Part Number: FLD9820-Clear

Dimensions: 98mmx20mm (2 steps)

Shade: Clear


NOBILCAM Flexible Disc

Technical Data




Flexural strength could denote the ability to resist deformation. NOBILCAM Flexible Discs provide higher strength than other flexible discs.

Flexural Strength




Elongation at Break



Elongation at break expresses the capability to resist changes of shape without crack formation. NOBILCAM Flexible Discs have better performance on softness and elasticity.

Excessive water absorption would affect the quality of flexible framework. According to the test results, NOBILCAM Flexible discs perform well.

Water Sorption



Form FLD Rev 11/2022