



## FULFILLING YOUR MILLING

### High Quality Zirconia Milling Discs



### GM-3D Gradient Multilayer

- Seamless gradient in translucency, strength and shade
- Nature Look, Nature Smiling
- Time saving and easy processing

Flexural strength	700-1050Mpa
Translucency	43%-57%
Sintered density	6.07±0.03g/cm <sup>3</sup>

### MHT Multilayer High Translucent

- All-in-one & One-for-all
- Seamless gradient makes colors change smoothly
- Time saving and easy processing

Flexural strength	1000 Mpa
Translucency	45%
Sintered density	6.07±0.03g/cm <sup>3</sup>

### HT High Translucent

- Exceptional strength with appealing translucency
- Wide indications for full contour
- Reproduce Vita\* shades perfectly
- Efficient and ideal results

\* VITA is a trademark of VITA ZAHNFABRIK H. RAUTER GMBH & CO.

Flexural strength	1200 Mpa
Translucency	43%
Sintered density	6.07±0.03g/cm <sup>3</sup>

# FULFILLING YOUR MILLING

## High Quality Zirconia Milling Discs

- ☉ We stock our milling discs in popular shades, sizes and thicknesses, shown below
- ☉ Lead time to order non-inventoried items is 2-6 weeks
- ☉ We may broaden our stock according to customer demand

### GM-3D Gradient Multilayer

GM	Part No. Prefix	Dimension	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D1	D2	D3	D4	BL1	BL2	BL3
	98GM3D14	98×14 mm																				
98GM3D16	98×16 mm																					
98GM3D18	98×18 mm																					
98GM3D20	98×20 mm																					
98GM3D25	98×25 mm																					

stocked item 

### MHT Multilayer High Translucency

MHT	Part No. Prefix	Dimension	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D1	D2	D3	D4
	98MHT14	98×14 mm																	
98MHT16	98×16 mm																		
98MHT18	98×18 mm																		
98MHT20	98×20 mm																		
98MHT25	98×25 mm																		

stocked item 

### HT High Translucent


HT	Part No. Prefix	Dimension	White	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D1	D2	D3	D4
	98HT12	98×12 mm																		
98HT14	98×14 mm																			
98HT16	98×16 mm																			
98HT18	98×18 mm																			
98HT20	98×20 mm																			

stocked item 

Part number= \*Part Number Prefix\* follow with Shade,e.g.98HT14-A3

Other sizes are available,please inquire.

Also available in following specifications per request



10mm - 25mm

**NOBILCAM™**