

All PMMA discs' flexural strength is over 120MPa

## **CAD/CAM Discs and Blocks**



### **Multilayer PMMA Discs and Blocks**

- ⦿ 20 shades available in stock
- ⦿ Natural tooth-like layering
- ⦿ Excellent wear resistance, shade stability and esthetics

### **Clear PMMA Discs and Blocks**

- ⦿ Crystal clear surface with outstanding translucency
- ⦿ Low water absorption and resists stains
- ⦿ Burns out completely without a trace



### **Monolayer PMMA Discs and Blocks**

- ⦿ 19 shades available in stock including clear and pink
- ⦿ Excellent wear resistance, shade stability and esthetics
- ⦿ Excellent finishing, bright surface

### **Full Denture PMMA Discs and Blocks**

- ⦿ Both single layer and multi layer are available in arch part
- ⦿ 4 pink shades available in base part
- ⦿ One-Step approach to full dentures



# CAD/CAM Discs and Blocks

## Product Availability

- © We stock our PMMA discs in popular shades, sizes and thicknesses, shown below
- © Lead time to order non inventoried items is 2-4 weeks
- © We may broaden our stock according to customer demand



### Monolayered PMMA Discs



Disc	Part No. prefix	SPEC.(mm)	A0	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	(0)Clear	2S	BL1	BL2	BL3	BL4	
	PMSL9810	98x10 0 Step																								
	PMSL9812	98x12 2 Steps																								
	PMSL9814	98x14 2 Steps																								
	PMSL9816	98x16 2 Steps																								
	PMSL9818	98x18 2 Steps																								
	PMSL9820	98x20 2 Steps																								
	PMSL9825	98x25 2 Steps																								
	PMSL9830	98x30 2 Steps																								
	PMSL9514	95x14																								
	PMSL9516	95x16																								
	PMSL9518	95x18																								
	PMSL9520	95x20																								
	PMSL9525	95x25																								

### Multilayered PMMA Discs

Disc	Part No. prefix	SPEC.(mm)	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	BL1	BL2	BL3	BL4				
	PMML9810	98x10 0 Step																								
	PMML9812	98x12 2 Steps																								
	PMML9814	98x14 2 Steps																								
	PMML9816	98x16 2 Steps																								
	PMML9818	98x18 2 Steps																								
	PMML9820	98x20 2 Steps																								
	PMML9825	98x25 2 Steps																								
	PMML9830	98x30 2 Steps																								
	PMML9514	95x14																								
	PMML9516	95x16																								
	PMML9518	95x18																								
	PMML9520	95x20																								
	PMML9525	95x25																								
	PMML9530	95x30																								

Stocking inventory available for quick shipment is indicated in pink  
 Part Number is combination of "Part No. Prefix" + "--Shade". Example: PMML9814-A1  
 10mm and 30mm thickness also available

Also available in the following specifications per request



Form PMMA Rev 09/2022