



### Natural Beauty

Life-like blend of shade and translucency  
 21 popular molds satisfy lab's needs  
 Familiar forms: Square and Tapering  
 Available popular shades: A1,A2,A3,B1

### High-performance property


Excellent shade stability  
 Low affinity to plaque

### Easy to Use


Easy to set-up posterior  
 No-wax carding, no need to wipe-off wax



### Physical Property



### Bonding Strength/N



### High Quality Raw Material

High quality raw material in use to achieve a reliable quality of the finished teeth products from the very beginning

### Reliable Shade Control System

Strict human verification  
 Second machine inspection



### Easy to Order

Website grid, phone, fax  
 Website bar code scanner linked to website  
 Digital order processing and fulfillment system





# NOBILDENT Select XL

Let Your Smile Shine with NOBILDENT



All measurements in millimeters



**SUGGESTED ARTICULATIONS**

UPPER ANTERIORS	LOWER ANTERIORS	POSTERIOR 28°/33°
T2	L2	30
T6	L6	32
T7	L7	32
T8	L8	34
T9	L9	*36M
T10	L10	*36M
S8	L8	34

**UPPER ANTERIOR TAPERING**



**UPPER ANTERIOR SQUARE**



**LOWER ANTERIORS**



**28° POSTERIOR**



**33° POSTERIOR**

