


A simple solution for non-powered therapy



- Sloped, soft, multi-density foam cradles heels to reduce pressure and shearing in this critical zone.

- Each air and foam filled cell communicates to create exceptional envelopment and pressure redistribution.
- Firm perimeter support for ease of transfer and safety



FUSION[®] G[™] Features

Dynamic Dispersion

Differs from all other non-powered alternatives. Each air and foam filled cell communicates to create exceptional envelopment and pressure redistribution.

Long Lasting Durability

Since the air is supporting a majority of the load, as opposed to the foam, the system is designed to last years beyond the competition.

Proven Redistribution Therapy

The "closed" system envelopes and redistributes pressure similar to other cellular based technologies. The competitor's "open" Self-Adjusting Technology vents making the single foam pieces handle the load individually.

Weight Limit: 300 lbs. (136 kg)

Mattress Size: 80" (or 84") x 36" 42" 48"x 6"


Cover: Water resistant, vapor-permeable, flame resistant and biocompatible

Warranty: One Year Limited Warranty on control unit and soft goods

LIT-0015 | R100121


Dynamic Dispersion™

how it works...




FUSION C™

Although, patient appears larger, points of contact are greater thus increasing pressure redistribution.



Outstanding envelopment and redistribution of pressure for a non-powered system

Competitive product shows higher peak pressures



The leading "Self-Adjusting Technology" competitor

