

The materials and information provided via this brochure are distributed internationally, but not all the products or services found on this brochure are necessarily available in your country or area. Please contact local office or agents to make sure of the availability of the relevant products or services.

**Mindray Medical Netherlands B.V.**  
Drs.Wvan Royenstraat 8, 3871 AN Hoevelaken, The Netherlands  
Voice: (31-33) 254-4911  
Fax: (31-33) 253-4280

**Mindray Medical Italy S.r.l.**  
Via Leonardo DA Vinci, 158, 20090 Trezzano Sul Naviglio (MI), Italy  
Voice: (39-02) 5737-401  
Fax: (39-02) 5523-2018

**Mindray Medical España S.L.**  
Avenida Manoteras 38, Edificio B, Bajo B008. 28050-Madrid Spain  
Voice: (34-91) 392-3754  
Fax: (34-91) 392-3756

**Mindray Medical Germany GmbH**  
Goebelstrasse 21, 3. OG, 64293 Darmstadt  
Voice: (49 6151) 3910 0  
Fax: (49 6151) 3910 300

**Mindray (UK) Limited**  
Mindray House, Kingfisher Way, Hinchingsbrooke Business Park, Huntingdon, Cambs, PE29 6FN  
Voice: (44) 1480-416-840  
Fax: (44) 1480-436-588

**Mindray Medical France SARL**  
Europarc Creteil 1 allée des cerisiers 94000 CRETEIL, France  
Voice: (33-1) 4513-9150  
Fax: (33-1) 4513-9151

**Mindray Medical Rus Co.Ltd**  
2 Zvenigorodskaya Street, 13, bld 41, 123022 Moscow, Russia  
Voice: (7-499) 553-6036  
Fax: (7-499) 553-6039

**Mindray Medical Colombia S.A.S**  
Room No.702, Torre UNIKA,Carrera 9 No.77-67, Bogota, Colombia  
Voice: (57-1) 3130892 / 3210916

**Mindray do Brasil – Comércio e Distribuição de Equipamentos Médicos Ltda.**  
Rua Tavares Bastos, 329, Perdizes, São Paulo - SP, CEP 05012-020  
Voice: (55-11) 3124-8026  
Fax: (55-11) 3078-8035

**Mindray Medical Mexico S. de R.L. de C.V.**  
Félix Parra # 175 Colonia. San José Insurgentes Sur Delegación Benito Juárez 03900 México, D. F.  
Voice: (52-55) 5661-9450 / (52-55) 5662-6620  
Fax: (52-55) 5662-6597

**PT.Mindray Medical Indonesia**  
The East 11th Floor, Unit 08 Jl. Dr. Ide Anak Agung Gde Agung Kav. E 3.2 No.1 Jakarta 12950  
Voice: (6221) 2902-7280  
Fax: (6221) 2902-7283

**Mindray Medical Egypt**  
Ground Floor -11 El Imam Metwally El Sharawi Street, Sheraton Heliopolis - El Nozha, Cairo, Egypt  
Voice: (2) 02 22691047  
Fax: (2) 02 22691046

**Mindray Medical Australia**  
ADD: 19 Harker St, Burwood VIC 3125, Australia  
Voice: (61-3) 9830 6610  
Fax: (61-2) 9550 2841

**Mindray Medical Nigeria**  
ADD: 1004 Estate,Victoria Island, Lagos, Nigeria  
Voice: (+234) 7030682053

**Mindray Medical India Pvt. Ltd. (Delhi Office)**  
16th floor, Building No. 9-B, DLF Cyber City, DLF Phase- 3, Gurgaon- 122002, Haryana  
Voice: (91-124) 4632-488  
Fax: (91-124) 4632499

**Mindray Medical India Pvt. Ltd. (Mumbai office)**  
B-404, City Point, Andheri Kurla Road, Andheri East, Mumbai, 400 059, India  
Voice: (91-22) 4020-0000  
Fax: (91-22) 4020-0011


**Mindray Medikal Teknoloji Istanbul Ltd. Sti.**  
KÖYALTI MEVKİLİ 29 EKİM CAD.N.1 KUYUMÇUKENT 2.NO.PLZ.3.KT. D.1- 9 34520 BAĞÇELİEVLER/YENİBOSNA, İSTANBUL, TÜRKİYE  
Voice: (90-212) 482 08 77  
Fax: (90-212) 482 08 78

**Mindray Medikal Teknoloji Istanbul Ltd. Sti. Branch Office Ankara**  
TURAN GÜNES BULVARI ALEKSANDER DÜBÇEK CAD. NO:18/A ÇANKAYA/ANKARA, TÜRKİYE  
Voice: (90-312) 441 96 92  
Fax: (90-312) 441 96 93

**Mindray Medical Vietnam Co.,Ltd.**  
Unit 2, 4th Floor, "B" Building, Broadway Office Park, 102 Nguyen Luong Bang Blvd., Tan Phu Ward, Dist. 7, HCMC, Vietnam  
Voice: (84) 8 5413-6970  
Fax: (84) 8 5413-6971

**Mindray Medical Thailand Limited**  
1768 Thai summit Tower Room 407, 4th floor, New petchaburee Road, Kwaeng Bangkapi, Khet Huai Khwang, Bangkok 10320, Thailand  
Voice: (66) 2652-7288  
Fax: (66) 2652-7266

**Mindray Medical Korea**  
KNI Tower 9F 823-19 Yeoksam dong, Seoul-ganum, Korea  
Tel: 008225688040  
Fax: 008225688042



# BeneHeart D6

Defibrillator/Monitor

Powerful life saving solution

Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China  
Tel: +86 755 8188 8998 Fax: +86 755 26582680  
E-mail: intl-market@mindray.com www.mindray.com

**mindray** | healthcare within reach are registered trademarks or trademarks owned by Shenzhen Mindray Bio-medical Electronics Co., LTD.  
© 2015 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved. Specifications subject to changes without prior notice.  
P/N:ENG-BeneHeart D6-210285x12P-20170727

**mindray**

**mindray**  
healthcare within reach




4-in-1 design and powerful function for full spectrum applications

With a 4-in-1 integrated design (manual defibrillation, AED, pacing, and monitoring modes), BeneHeart D6 puts any unexpected circumstances under your control.


#### Manual Defibrillation

Asynchronised defibrillation mode for cardioversion of ventricular fibrillation. Synchronised defibrillation mode for cardioversion of atrial fibrillation.



#### Multi-parameter Monitoring

BeneHeart D6 offers comprehensive monitoring functions designed to give clinicians the information they need. In addition to ECG with respiration, clinicians can choose 12-lead ECG with full reports, SpO<sub>2</sub>, 2xTemp, NIBP, 2xIBP and EtCO<sub>2</sub>.



#### AED

In AED mode, BeneHeart D6 automatically analyses the rhythm and indicates whether or not a shockable rhythm is detected. Voice and text prompts guide the user through the process. Voice recording (180 minutes) is also available for case review.



#### Non-invasive Pacing

BeneHeart D6 offers demand mode and fixed mode pacing with adjustable rates and output. The 4:1 key enables clinicians to quickly select 1/4 of the pacer rate for observation of the underlying rhythm.




# Fast defibrillation

## The fastest defibrillator

Mindray strives for constant innovation to improve the clinical aspects of product performance. The new generation of technology platform enables Mindray to improve the performance of the BeneHeart D6 defibrillator to meet changing clinician needs.

BeneHeart D6 gives you a greater chance of success for those patients suffering cardiac arrest. It only takes 7.5 seconds to complete the whole defibrillation operation. Studies show that when a patient suffers cardiac arrest, success rates for defibrillation drop for every second between CPR and defibrillation shock. Every second counts for cardiac arrest patients.\*



Manual defibrillation with clear 1-2-3 steps



1. Select Energy
2. Charge
3. Shock



External Paddles with function buttons

Buttons for energy selection, charging and shock delivery improve usability for clinicians.


\*Edelson DP, Abella BS, Kramer-Johansen J, et al. Effects of compression depth and pre-shock pauses predict defibrillation failure during cardiac arrest. Resuscitation. 2006 Nov;71(2):137-45.

# Effective resuscitation



## 360J high energy

BeneHeart D6 defibrillator/monitor features 360J biphasic technology, which increases the chance to save difficult-to-defibrillate patients. Studies have shown that cardiac arrest is common among ventricular fibrillation (VF) patients and that defibrillation of recurring episodes of VF is increasingly difficult. A randomised controlled clinical trial shows the rate of VF termination increases with charge energy, when charge energy is 200J and above.\*



When impedance is adjusted, voltage and energy delivery cycle are adjusted automatically to correspond with impedance

\*Stiell I, Walker R, Nesbitt L, et al. Biphasic Trial: A randomized comparison of fixed lower versus escalating higher energy levels for defibrillation in out-of-hospital cardiac arrest. *Circulation*. 2007;115:1511-1517.


## CPR solution

### Accordance with the 2015 AHA/ERC guideline

The CPR Sensor is easy to use with its ergonomic design and graphic mark of the pressure spot. It delivers instant audiovisual feedback so that compressions are of proper depth and rate. It indicates complete chest recoil, hands-off time to help improve the patient's chance to survive.

## Intuitive contact impedance indicator


Colour coded indicator with real contact impedance value provides a more intuitive guide to clinicians.



CCI (CPR compression index): display dynamic change of CPR compression quality, including information such as compression rate, compression depth, compression recoil, compression fraction


Compression waveform area

Compression rate area & Compression depth area




## Adult/Paediatric mode

When changing from adult to paediatric mode, the default shock energy, monitoring range and parameter alarm limits change automatically to deliver the best treatment effect for both types of patients.



# Professional pre-hospital solution

## Pre-Hospital Emergency Information System(PHEIS)



### Tracking the Entire Resuscitation Process

PHEIS sends real-time patient data, such as ECG, SpO<sub>2</sub>, EtCO<sub>2</sub>, 2 Temp and NIBP to the hospital, allowing doctors to track the entire rescue process all the way from the field to the hospital.

### Remote Diagnosis to Provide Rapid Patient Treatment

For patients in a critical condition, every second counts. PHEIS helps to save valuable time by supporting remote diagnosis and allowing the hospital to prepare treatment well in advance. The Glasgow 12-lead ECG report helps doctors quickly decide on the reperfusion method (intervention or fibrinolytic therapy) for a suspected STEMI (ST-elevation myocardial infarction) patient.



## Reliable in any emergency environment

BeneHeart D6 can be conveniently charged inside the vehicle. While the optional secondary battery doubles its running time, the external charger carried in the vehicle can be used to charge any additional backup batteries. BeneHeart D6 complies with IP44 dust and water proof standard, which allows it to work in a variety of complex outdoor environments. The carrying case with cover shield protects D6 from knocks and the expansion bags can be used to store extra accessories.



High performance battery & Battery charger

Excellent battery performance guarantees you can deliver defibrillation at any time. LED indicator shows battery status clearly. Battery charger keeps two spare batteries in charge and ready to use.




Vehicle Docking station

Docking station with built-in power supply simplifies the anchoring and locking process of the defibrillator.



Carrying case


The large capacity carrying case adds protection to the BeneHeart D6 and eases the transport of accessories.



## Effective IT solution

Our simple yet effective IT solution manages all the information from BeneHeart D3 defibrillator/monitors to avoid manual recording and so improve efficiency and reduce the workload of clinical staff. All information can be databased.

A simple yet robust network connection following standard information infrastructure in most hospitals: Transmit data through 5G/2.4G WiFi, international standard IHE HL7 protocol and DHCP to obtain IP address automatically.




Monitor information transmits to CMS

Centralised management of monitoring information improves efficiency.



Rescue information transmits to CMS

Patients treatment process can be reviewed conveniently at any time, providing data support for further study.



Device status transmits to CMS  
All the testing results including self test and user test will be sent to Mindray CMS. Centralised device management and maintenance ensures every defibrillator is "Ready" to use.