

## WHICH SIZE TIDY TUBING IS RIGHT FOR YOU?

	1833	1834	1835
Extended Length	5 Feet	10 Feet	15 Feet
Recommended Use	<ul> <li>Portable O2 Concentrator</li> <li>Oxygen Tank Bag</li> </ul>	<ul> <li>Portable O2 Concentrator</li> <li>Pull Cart</li> <li>O2 Tank Bag</li> <li>Wheelchair</li> </ul>	<ul> <li>Stationary Concentrator</li> <li>Wheelchair</li> <li>Pull cart</li> <li>Portable Concentrator (when stationary)</li> </ul>
Included	<ul> <li>Clothing Clip (prevents nose tugging)</li> <li>2 foot Vanish Cannula</li> <li>Swivel connector (prevents tangling)</li> </ul>	<ul> <li>Clothing Clip (prevents</li> <li>nose tugging)</li> <li>2 foot Vanish Cannula</li> <li>Swivel connector(prevents tangling)</li> </ul>	<ul> <li>Clothing Clip (prevents nose tugging)</li> <li>2 foot Vanish Cannula</li> <li>Swivel connector (prevents tangling)</li> </ul>

- Connect up to three 1835's for a recommended maximum of 45'
- Annual replacement recommended
- Install clothing clip near waist to prevent nose tugging
- Recommended use will vary by person and oxygen supply system

## www.captivetechnologies.com

