



Cost effective care without compromise. The Joerns WeCare™ bed continues our legacy of quality and reliability, now in a convenient package for easy delivery to both institutions and community. Designed and field tested for one-person delivery and setup, WeCare brings fall prevention, resident comfort and caregiver safety to any setting.



Enhanced Capabilities and New Options

- **UNPARALLELED STABILITY:** Enhance confidence entering / exiting the bed with full range stability and foot placement.
- **A SAFER SYSTEM:** Protect against falls injury and reduce caregiver strain with the 8.0" ultra-low height and ergonomic 28" high height.* Add peace of mind knowing WeCare is guaranteed to meet bed entrapment guidelines.
- **ENHANCED COMFORT AND CARE:** Augment care, improve sleep, and reduce movement in bed with a fully adjustable sleep deck, AutoContur positioning, and optional staff control with Trendelenburg/Reverse Trendelenburg and function lock.
- **A BETTER FIT:** Improve comfort, simplify repositioning and reduce falls risk with flexible sizing. WeCare offers optional 4" length extension and width extensions from 35" to 39" (Twin XL) and 42".
- **EASIER INSTALLATION:** Protect installers from excessive strain and simplify delivery with WeCare's unique transport dolly and stair climber.
- **FLEXIBLE DELIVERY:** Receive WeCare the way you want it. Packaged for easy delivery to suit your needs, the bed system can ship palletted or by parcel carrier.



Standard Offering

- ◊ **Reduced risk of injury:** With greater travel range
- ◊ **Exceptional clearance:** Open architecture legs
- ◊ **Comfort and shear reduction:** Three-section mattress support platform
- ◊ **Better sleep deck:** Therapeutic grid surface
- ◊ **Industry leading stability:** Durable tubular steel design

- ◊ **Intuitive ergonomic hand pendant:** 8-function control with integrated hook
- ◊ **Safer repositioning and egress:** Entrapment compliant assist device can be mounted at head and foot, plus optional staff control offers gravity assisted repositioning
- ◊ **Effortless operation:** Four quality DC motors
- ◊ **Managed controls:** Pendant function lock-out available with optional staff control

*Basic Bed: 8.75" - 26"

Delivering care together...



We Care about delivering fast and effective care as much as you do! We've made ordering and setup of your bed quick and easy. Receive WeCare the way you want it. Packaged for easy delivery to suit your needs, with flexible delivery options. Joerns WeCare elevates patient comfort and care with simplified care solutions.

Extending a helping hand



F028 Deluxe Assist Handle



F025 Assist Handle



F026 Assist Handle



F14 Assist Device



F168 Assist Device
(requires N832AL mount)

Caring for more patients



Length Extension



Width Extension



Multiple Panel Style
and Finish Options



Trapeze Patient
Helper

Delivering care safely



Transport Dolly
and Stair Climber



Storage and
Transport Cart

Ordering Part Numbers

WECC-WECARE-RIL	WECC-WECARE-RAH	WECC-WECARE-BASE
Roll-in-Low Bed Split-frame Design 8" – 28" Height 80" Length	Roll-Any-Height Bed Split-frame Design 8" – 28" Height 80" Length	Base Bed Split-frame Design 8.75" – 26" Height 80" Length

Specifications

Bed Mobility	Roll-in-Low / Roll-any-Height	Casters	Four casters (two locking)
Hand Pendant Functions	Standard: 3 Function	Safe Working Load	450 lb 400 lb Max Patient Weight
	Plus AutoContour	Underbed Clearance	18"
Low Height	8.0" (Base Model 8.75")	Product Weight	145 lb (77.7 lb foot / 67.3 head)
High Height	28" (Base Model 26")	Electrical	120 VAC, 60 Hz, 4.0 Amps, Class II
Bed Length	Standard: 80"	Regulation Compliance	Conforms to IEC 60601-2-52
	Optional: 4" Extender	Warranty	1 year mechanical / electrical, 7 year structure, life welds
Bed Width	Standard: 35"		
	Optional: CareWide Accessory 39" (TwinXL) and 42"		