Call us at 1-800-944-1997 | sales@evrider.com | www.evrider.com

There is Power Mobility. And then there's

POWERFUL MOBILITY



There is Power Mobility. And then there's

POWERFUL **MOBILITY**

Move-X Rollator

The most compact Rollator in the world. This latest generation of Rollator is amazingly light, comfortable, and durable. The Move-X folding in a second on your index fingers makes the Rollator extremely compact for transportation. This work of art is like a luxury car, which makes you feel confident using it. With Move-X, you will no longer find it difficult to go out. It's a beautiful day, go out and enjoy. Move to a wonderful life.







LOAD CAPACITY



BATTERY



RANGE



Handle Width	19 in / 47 cm
Overall Height	30.5 in / 77.5 cm to 35.5 in / 89.5 cm
Folded Overall	25 in x 11 in x 15 in 64 cm x 28 cm x 39 cm



Folded Overall Width	24 in / 61 cm
Weight	16 lbs / 7.6 kg
Seat Height	22.8 in / 58 cm



Seat Width	15 in / 38 cm
Outside Diameter	8 in / 20.3 cm
Load Capacity	250 lbs / 113.4 kg
Colors	Black, Blue