

# Specifications

Colour	Model Numbers	
<span style="color: #4F81BD;">■</span> Sky Blue <span style="color: #002060;">■</span> Midnight Blue	600-210006 600-210007	600-190006 600-190007
<b>Overall width</b>	24.75" / 63 cm	24.75" / 63 cm
<b>Folded width</b>	10.5" / 26.7 cm	10.5" / 26.7 cm
<b>Length</b>	28.25" / 71.75 cm	28.25" / 71.75 cm
<b>Handle height range</b>	32.4" – 40.2" 82.3 cm – 102 cm	29.4" – 37.2" 74.7 cm – 94.5 cm
<b>Distance between handles</b>	17.75" / 45 cm	17.75" / 45 cm
<b>Seat Height</b>	21" / 54 cm	19" / 48 cm
<b>Seat width</b>	16" / 40.6 cm	16" / 40.6 cm
<b>Seat depth (rollator)</b>	12.5" / 31.75 cm	12.5" / 31.75 cm
<b>Seat depth (transport chair)</b>	16" / 40.6 cm	16" / 40.6 cm
<b>Weight (Transport)</b>	30 lbs (13.6 kg)	30 lbs (13.6 kg)
<b>Weight (Total)</b>	32 lbs (14.5 kg)	32 lbs (14.5 kg)
<b>Weight Capacity</b>	275 lbs (125 kg)	275 lbs (125 kg)



SKY BLUE



MIDNIGHT BLUE




Limited lifetime warranty

Please visit:

[triumphmobility.com](http://triumphmobility.com)

for additional information



SCAN ME!

## Contact us

Triumph Mobility

151 Bentley Ave, Unit 4

Nepean, Ontario, K2E 6T7

Tel: +1 (855) 546 0711

Fax: +1 (866) 225 3059

Email: [info@triumphmobility.com](mailto:info@triumphmobility.com)

[www.triumphmobility.com](http://www.triumphmobility.com)



TriumphMobility

©2020Triumph Mobility Inc. All rights reserved.  
 Triumph™ and Prestige™ are trademarks of  
 Triumph Mobility Inc. The Triumph Prestige  
 rollator and transport chair is covered by pending  
 patent and industrial design protection.

Printed in Canada

File name: R\_202204\_TriumphPrestige\_brochure\_EN

Proudly designed in Canada



**TRIUMPH**  
**prestige**

# The Prestige Lifestyle

- ✓ No tools needed
- ✓ No fussy instructions

Finally, a support to keep pace with your busy life.



ROLLATOR MODE

TRANSPORT CHAIR MODE



## Perfect Flexibility

10.5 in



The Prestige locks closed at **10.5 in** at its widest point – perfect for tucking it away in a restaurant or at home. Its economic size also makes it easy to fit in a car trunk, weighing only **13.6 kg**



The Prestige transforms from rollator to transport chair in seconds: **swivel the flexible padded back rest, flip the handgrips around**, you're ready to move!



Don't worry about wearing it out: the frame has a **lifetime warranty** and all other parts are guaranteed for 2 years.



≤ 11 lbs (5 kg)

The basket underneath makes the Prestige perfect for shopping, securely attached at six different points and strong enough to **hold up to 11 lbs.**

## What else?

### Comfortable



- The padded backrest adjusts to the natural angle in both modes
- A padded seat cushion with a roomy seat depth of 12.5" in rollator mode and 16" in transport chair mode
- Ergonomic, angled and adjustable handgrips
- Padded armrest at a natural position of 9.5" from seat to arm

### Stylish



- Two fashionable colours to choose from
- Wires all internal for a sleek look with nothing to snag

### Safe



- Soft tires provide a natural suspension and balance
- Drum brakes provide reliable all-weather braking
- Sturdy aluminum frame design