

FreeRider[®]USA



Smart Fold, Easy Go

Luggie Classic

LENGTH	39 in (99.1 cm)
WIDTH	18 in (45.7 cm)
LENGTH (FOLDED)	26 in (66 cm)
WEIGHT CAPACITY	250 lbs. (113.4 kg)
WEIGHT (SCOOTER)	51.7 lbs (23.5 kg)
WEIGHT (SCOOTER+BATTERIES)	71.4 lbs (32.4 kg)
BATTERY	2x 12AH Lead Acid
BATTERY CHARGER	2-amp off-board charger
MAX. SPEED	4 mph (6.4 kph)
RANGE*	10 miles (16.1 km)
TURNING RADIUS	35.5 in (90.2 cm)
GROUND CLEARANCE	2.5 in (6.4 cm)
MAX. INCLINE	6 degrees
WHEELS/TIRES	Front: 2" x 6" PU solid Rear: 2" x 7" PU solid
DRIVE SYSTEM	24-volt DC motor, sealed mini-transaxle, rear-wheel drive

STANDARD COLORS ● ●

The information contained herein is correct at the time of publication. Factors such as battery charge, battery condition, user weight, tire pressure, and terrain type will cause variation to maximum range and speed. Due to continual product improvement and manufacturing tolerances; length, width, turning radius and ground clearance can be subject to variance of +/- 3%. We reserve the right to alter specifications without prior notice.



Economical • Compact • Simple

Easy-folding mobility scooter with tight turning radius for effortless maneuverability. Conveniently designed for daily use.



Toll Free [888] 957-1888 | www.FreeRiderUSA.com