




M07340-1




STEP 1:
Attach front part to rear part




Brake cable




STEP 2:
Attach wheel to front part



STEP 3:
Fix brake cable with front wheel




STEP 4:
Attach basket to front part




REV
00

1.	2.	3.	4.	5.	6.	7.	8.	9.	10.	11.
1x	1x	1x	1x Ø 310	1x M12x200	2x Ø 12	1x M12	2x M6x15	2x Ø 6	2x M6	1x, 1x, 2x


M07340-1




STEP 1:
Attach front part to rear part




Brake cable




STEP 2:
Attach wheel to front part



STEP 3:
Fix brake cable with front wheel



STEP 4:
Attach basket to front part



1.	2.	3.	4.	5.	6.	7.	8.	9.	10.	11.
1x	1x	1x	1x Ø 310	1x M12x200	2x Ø 12	1x M12	2x M6x15	2x Ø 6	2x M6	1x, 1x, 2x