


## Pessary Selection Chart

Symptoms	Prolapse 1st Deg.	Prolapse 2nd Deg.	Prolapse 3rd Deg.	Cystocele	Rectocele	Incontinence	Displaced (Retro)
Pessary Type							
RING	X	X		X			
CUBE		X	X	X	X		
GELLHORN		X	X	X			
SHORT-STEM GELLHORN		X	X	X			
DONUT		X	X	X	X		
SHAATZ	X	X	X				
GEHRUNG		X	X	X	X		
DISH	X	X	X	X		X	
CUP	X	X	X	X		X	
OVAL	X	X		X			
RING WITH KNOB	X	X		X		X	
MARLAND	X	X	X	X		X	
HODGE	X	X		X		X	X
INCONTINENCE RING						X	
FLEXISHELF		X	X				

**\*Clinics, Pharmacies and other healthcare establishments can create account online and order at [www.canadapessary.ca](http://www.canadapessary.ca)**

Disclaimer: This page is for information purpose only and should not be used as a professional resource or advice to diagnose, treat or manage any health conditions. The information, including but not limited to text, graphics, images and other material contained on this page are for information purposes only and do not replace instructions provided by a Healthcare provider. No material on this site is intended to be a substitute for professional medical advice, diagnosis or treatment.