

DeVilbiss 5 Liter Oxygen Concentrator

HCPCS: E1390



FIGURE A

- Built upon a legacy of reliable, field-proven concentrator models, the DeVilbiss 5 Liter Oxygen Concentrator was designed with enhanced quality, durability and simplicity in mind
- Patented DeVilbiss Turn-Down Technology minimizes wear on internal components, reduces power consumption by 15% or more and extends the life of the concentrator
- Exclusive DeVilbiss OSD (Oxygen Sensing Device), standard on every unit, ensures patient safety and reliability for longer service intervals
- Readily accessible patient controls, protected cannula fitting and recessed humidifier nook to prevent damage
- Alarms Audible and visual high/low pressure, low flow, low oxygen, power fail, Oxygen Sensing Device
- Sold separately, the Transfill Caddy is designed with a steel frame to simply connect the 18 x 14 inch shelf, designed to support up to 40 lbs, to the DeVilbiss 5L Concentrator. Rear view (Figure A), Front View (Figure B)

PRODUCTS

Description

525DS	5 Liter Compact Concentrator
525DS-EW	525DS Concentrator w/5yr Warranty, 1/ea
525DS-Q	5 Liter Compact Concentrator Quiet
525DS-QEW	525DS CncntrtrQuiet w/5yr Warranty, 1/e

SPECIFICATIONS

Description

Delivery Rate	0.5 to 5 LPM
Outlet Pressure	8.5 psig
Oxygen Percentage	1-5 LPM, 87% - 96%
Power Consumption	310 Watt Average, 275 Watts @ 1.2 LPM and Below
Product Dimensions	24.5" (H) x 13.5" (W) x 12" (D)
Product Weight	36 lbs
Warranty	3 Years

ACCESSORIES

Description

303DZ-605	AIR FILTER 6-PKG - DFT, PKG
MC44D-605	EXT LIFE INTAKE BACT FILTR-DFT, EA
PV5LD-651	FINAL BACTERIA FILTER - DFT, 1/EA, EA
515LF-607	FLOW METER PACKAGE, LOW OUTPUT, EA
525DD-650	Transfill Rack System, 525DS Series, 1ea, EA
PF1100TUB	Transfill Connection Tubing, 1/ea, EA

PARTS

Description

525DD-626	Comp Filter Assbly Sintered, EA
-----------	---------------------------------

PARTS	Description
525DD-635	Aux Port SS Pkg, EA
525DD-603	Base Assembly, 1/ea, EA
525DD-636	Plastite Screw SS Pkg, EA
525D-613	CIRCUIT BREAKER, EA
525D-625	COMPRESSOR S.S.P., EA
525DD-625	Compressor box, EA
525D-643	COMPRESSOR REBUILD KIT, EA
515DZ-634	COOLING FAN PACKAGE, EA
515A-705	EXHAUST MUFFLER PKG, EA
525DD-639	DOOR INTAKE FILTER PACKAGE, EA
525DD-611	Front Cover Package, 1/ea, EA
525DD-627	Front label, EA
PV5LD-617	HOUR METER - DFT, EA
525D-615	LIGHT PANEL S.S.P., EA
505IZ-609	MOTOR MOUNT PKG 4-PKG - DFT, EA
525DD-606	Concentrator (5L, 10L) Outlet Port, 1/ea, EA
525DS-622	PC BOARD ASSY, EA
525DD-612	Rear Cover Package - No Aux Port, 1/ea, EA
515ADZ-702	ROTARY VALVE REPLACEMENT, EA