

# Benchtop Multifunction Shaker

## Horizontal + Pendulum + 3D



ID: PLP-DS-3D100

### Features

- 3 Independent modes: Swing, Horizontal and 3D
- In 3D mode, the mixing intensity can be set separately
- Automatically stops in horizontal position after operation
- Timer and power-off memory function
- Large LCD Screen for easy operation
- Equipped with a stepper motor for precise speed control and long lifespan
- Rust-proof design on the base plate for durability
- Ideal for a modern lab in need of a compact shaking or mixing unit



### Specifications

Speed Range	Horizontal Mode: 5-100rpm Pendulum Mode: 5-100rpm 3D Mode: 5-100rpm
Radius of Gyration	In Plane Mode: 15mm Swing Angle in Swing Mode: 10 degrees Speed Ratio 3D Mode: 10%-100%
Tray Size	338 x 308mm
Max Load	5kg
Timer	0-999h 59min
Power	35W
Voltage	AC 240V 50Hz/60Hz