

Hotplate Stirrer with Large Ceramic Plate

ID: PLP-SMS-200A



Features

- Heating and stirring with intelligent speeding function
- Stirring volume capacity up to 10 litres
- Temperature sensor PT100 with 440mm rod and clamp
- LED display for temperature and speed range
- Automatic stirring bar homing, magnetic stir bar area
- Anti-scalding design, flame-retardant
- Over temperature protection



Magnetic stir bar collection area



Automatic stirring bar homing



Specifications

Speed Range	100 - 2000rpm
Speed Control Accuracy	100rpm
Heating Temp. Range	RT-300 °C
Temp. Control Accuracy	Up to ±1 °C
Overheating Protection	320 °C
Max. Stirring Capacity	10L
Platform Size	150mm × 150mm
Platform Material	Cast Aluminum Sprayed Ceramics
Timing Range	1 - 999h 60min
Input Voltage	220V 50/60Hz
Net Weight	2.4kg