

Horizontal Shaker with Spring Clamp Platform

ID: PLP-SYC-WP200



Features

- Spring grid design, providing stable sample placement
- Dual-layer spring bracket allows for holding shorter vessels like culture dishes and taller vessels like conical flasks
- Automatic leveling to prevent sample spillage
- Power-off memory function, automatically resumes operation after power is restored.
- Stepper motor for precise speed control and long service life
- LED display for speed and time

Specifications

Speed Range	5-200rpm
Radius of Gyration	15mm
Max. Load	5kg
Timing Range	0-999 H 59 min
Power	20W
Voltage	AC 240V 50Hz/60Hz
Platform Size	380 x 330mm
Weight	9.6kg

