

Vacuum Aspirator Large Capacity


ID: PLP-SZY-2L



The vacuum aspirator is capable of collecting and retaining a range of liquids, including different types of cell culture dishes, test tube lids, and even experimental workstations, ensuring the laboratory remains safe and organized.

Features

- Compact benchtop instrument
- Automatic overflow switch, preventing overflow
- Two Liquid control methods:
Foot Pedal and knob switch
- Adjustable speed suction
- Quick connect fittings ensure secure connection to the various adapters



Specifications

Max. Suction Speed	1.6L/min
Air Filtration	0.45 µm Hydrophobic Filtration Membrane
Ambient Moisture	5-40°C 80%RH
Pipetting Bottle Volume	2L
Noise	≤60 dB
Power Display	LED Display
Battery Capacity	2600 mAh. Rechargeable
Power	6W
Power Adapter	Input 100-240V AC50/60Hz, 24V 1A
Product Dimension	147*147*286mm
Weight	0.8kg

