




## Microplate Centrifuge





## **Features**

- Can accommodate up to 2 microplates
- The motor has brakes to achieve a slow stop in 10s
- Inductive lid switch, which stops the motor once the lid is opened
- Bottom plate counterweight to prevent vibration







## **Specifications**

•	
Capacity	2 x 96-well PCR plate
Max. Relative Centrifugal Force	550 g
Motor	DC brushless
Rotor	Vertical fix
Rotating Speed	2,200 rpm
Noise	≤55 dB
Deceleration Time	≤10s (2,200 - 0 rpm)
Acceleration Time	≤10s (0 - 2,200 rpm)
Dimensions	256 × 224 × 203 mm
Voltage	AC 100 ~ 240V, 50 ~ 60HZ
Weight	3.6 kg