

AHS Disposable Nitrile Exam Gloves provide a safe and reliable solution for healthcare environments. The ambidextrous nitrile exam gloves are powder-free, chemo-rated and have a thickness of 3.5mm. Available in sizes from Small to Extra Large, the gloves come in boxes of 100 or cases 1000 (10 boxes). Keep your staff safe and protected with these tried-and-tested American Hospital Supply gloves.



Item Code	Description	Quantity	Barcode GTIN	Packaging Dimensions
AHS-GN-S_BX	AHS Nitrile Gloves Chemo-Rated Small	Box of 100 Each	850046814648	9" x 5" x 2.4"
AHS-GN-S_CS	AHS Nitrile Gloves Chemo-Rated Small	Case of 1000 Each	850046814655	12.4" x 10.1" x 9.6"
AHS-GN-M_BX	AHS Nitrile Gloves Chemo-Rated Medium	Box of 100 Each	850046814686	9" x 5" x 2.4"
AHS-GN-M_CS	AHS Nitrile Gloves Chemo-Rated Medium	Case of 1000 Each	850046814693	12.4" x 10.1" x 9.6"
AHS-GN-L_BX	AHS Nitrile Gloves Chemo-Rated Large	Box of 100 Each	850046814716	9" x 5" x 2.4"
AHS-GN-L_CS	AHS Nitrile Gloves Chemo-Rated Large	Case of 1000 Each	850046814723	12.4" x 10.1" x 9.6"
AHS-GN-XL_BX	AHS Nitrile Gloves Chemo-Rated XL-Large	Box of 100 Each	850046814747	9" x 5" x 2.4"
AHS-GN-XL_CS	AHS Nitrile Gloves Chemo-Rated XL-Large	Case of 1000 Each	850046814754	12.4" x 10.1" x 9.6"

- Powder Free
- Ambidextrous
- Medical Grade
- Chemotherapy Tested per ASTM D6978.
- Users should review MSDS for chemotherapy drugs prior to use to ensure level of protection required.
- Storage Recommendations: Avoid freezing and excessively hot temperatures. Keep away from moisture, florescent lights, x-rays, direct sunlight, and ozone.
- Meets ASTM D6319



Single Use Only.
Do not reuse



Keep Dry



Not Made with
Natural Rubber Latex



Date of Manufacture



Nonsterile



Lot #

CAUTION: BE SURE TO REVIEW DRUG LABELING AND SDS IF GLOVES ARE USED FOR PROTECTION AGAINST CHEMOTHERAPY DRUG EXPOSURE TO ENSURE THE GLOVES OFFER THE SPECIFIC TYPE OF CHEMICALS USED. REVIEW DRUG LABELING OR SAFETY DATA SHEETS FOR THE CHEMICALS BEING USED TO DETERMINE AN ADEQUATE LEVEL OF PROTECTION FOR THE INTENDED USE.