

# CLEANING & CARE GUIDELINES

INSTRUCTIONS OF USE

PHILLIPS-SAFETY.COM

COPYRIGHT ©2025. ALL RIGHTS RESERVED



**PHILLIPS SAFETY**  
PRODUCTS INC.

## USE

---

Phillips Safety Products must be use for it the designated propous only.

## WARRANTY

---

Phillips Safety Products have one year warranty against manufacturer defects. Attention: This does NOT include normal wear and overloading of use.

## CLEANING

---

Use a mild detergent or cleaner to remove any dirt or debris.  
Do NOT use abrasives.

## DISINFECTING

---


It can be disinfected with most standard hospital grade disinfectants containing: QATs, IPA, or Hydrogen Peroxide.

# CONTACT

Should you need any further information,  
please do not hesitate to contact us.

123 Lincoln Blvd, Middlesex, NJ 08845, USA 

+1 (866) 757 1307 

service@phillips-safety.com 

www.phillips-safety.com 