


SAY

Hello
TO...



BETTY &
BERTIE
MOBILITY

INFO@BETTYANDBERTIE.COM

0800 999 6606

UNIT G, THE COURTYARD BUSINESS CENTRE,
LONESOME LANE, REIGATE, SURREY, RH2 7QT

HAPPY CUSTOMERS

You're in good hands with us! Here are some of our top reviews by our customers which we think you'll enjoy...



There are no ends to where you can go with this scooter, it is very smart and compact meaning I can leave it in the boot of my car ready to go. I think the controls are very simple to use and the automatic folding is amazing!

- Carol & Buster



Asked all the questions I needed to ask, they took the time to understand my situation rather than pushing a sale on me and explain how my scooter would work. I was made to feel important and valuable with Betty and

Bertie

- Rose Jay



I love the colour of the black and orange combination, the brightness of them is very attractive. I think the opening with the touch of a button is brilliant and the controls are very efficient which helps my arthritis

- Pog & Jim



When your beautiful scooter arrived safely, well packed and in absolutely 1st class condition it was like re- charging my pacemaker! Thank you for a splendid product, 1st class service. I would certainly recommend Betty and

Bertie again.

- Patrick Mansen

Check out our...

**MANUAL
FOLDING**

MOBILITY SCOOTER

LUPIN



UK'S MOST CONVENIENT LIGHTWEIGHT FOLDING MOBILITY SCOOTER



Travel

Lupin is incredibly travel friendly. Weighing only 23kg and folding compactly, it can fit snugly into the back of your car and easily lifted in and out of taxis, trains and busses. The lithium battery is also airline safe so it can travel on holiday with you too.



Safety

Lupin comes with a rear-draw bar with anti-tip wheels, allowing you to store bags or luggage safely. Its strong, flat foot-bed makes getting on and off the Lupin easy and safe. Its powerful disc brake, LED headlight and wide rear axle ensure safety at all times.



Convenience

Lupin's lithium-ion battery can be removed and connected to a power outlet for hassle-free charging, and offers a huge travel range of up to 24km from a single charge. No assembly required; just safety, comfort and convenience all at your fingertips.

KEY BENEFITS



AIRLINE FRIENDLY



TRAVEL FRIENDLY



ONLY 23KG



EASY TO USE



SMART BRAKING SYSTEM



LIGHTWEIGHT & COMPACT

Simply pull a handle and lift the seat, and watch it unfold from being compact and transportable, to a comfortable and simple-to-use solution to all your portable mobility scooter needs!

EASY TO USE

1



2



3



SPECIFICATION

BODY

- Dimensions (Unfolded):
H: 87cm W: 56cm L: 98cm
- Dimensions (Folded):
H: 54 cm W: 56cm L: 98 cm
- Ground clearance: 4.7 inches / 12cm

WHEELS

- Wheel size: diameter F: 25 cm / R: 25 cm
- Wheel type: Pneumatic
- Maximum incline: 8 degrees +-1 degree

BATTERIES

- 1 removable lithium-ion battery that goes for 30 - 35km / 18m - 21m on a single charge
- Charging time: 10 - 12 hours
- Airline friendly (when used with airline friendly battery)

KEY FEATURES

- ✓ Weighs Only 23kg
- ✓ Load Capacity: 100kg / 15.7 stone
- ✓ Small Enough To Fit Into The Boot Of A Car
- ✓ Safety Light
- ✓ Anti Tip Wheels
- ✓ Adjustable Handlebar
- ✓ Intelligent Braking System
- ✓ Easy Speed Control
- ✓ Removable Battery
- ✓ Expandable Footrest
- ✓ Luggage Carrier
- ✓ Easy To Fold & Store

FAQ

How do I claim back VAT?

Purchase your scooter and give us a call or drop us an email and we can refund the VAT back to you if you can provide us with a reason for your VAT relief.

How does this scooter controlled?

This scooter is operated with a button control panel found on the top of the handlebar. Ensure the scooter is on and a mode is selected to use.

Does this scooter come in other colours?

Yes this scooter does come in other colours. We offer a Black/Red or a Grey/Black design of this model.

How do I charge the battery?

Remove the battery from the scooter and connect your charger to the port and plug the mains into the wall. Batteries will take around 10-12 hours to fully charge.

How can I put this scooter in my car?

Ensure the scooter is folded down and ensure to lift the scooter from the base as this will provide the most grip and distribute the weight evenly.

Is this scooter waterproof?

This scooter is not waterproof and we suggest drying it off with a cloth and removing the battery if in contact with heavy rain.

£1620 INC VAT

£1350 EX VAT

REN



WHAT IS THE DIFFERENCE BETWEEN REN & LUPIN?

We have made this model Ren much stronger as it has a higher **load capacity of 165kg (25st)** AND

Ren has **dual motors** making it so much better on hills and tougher terrain!

UK'S MOST CONVENIENT LIGHTWEIGHT FOLDING MOBILITY SCOOTER



Travel

Ren is incredibly travel friendly. Weighing only 23kg and folding compactly, it can fit snugly into the back of your car and easily lifted in and out of taxis, trains and busses. The lithium battery is also airline safe so it can travel on holiday with you too.



Safety

Ren comes with a rear-draw bar with anti-tip wheels, allowing you to store bags or luggage safely. Its strong, flat foot-bed makes getting on and off the Ren easy and safe. Its powerful disc brakes, LED headlight and wide rear axle ensure safety at all times.



Convenience

Ren lithium-ion battery can be removed and connected to a power outlet for hassle-free charging, and offers a huge travel range of up to 24km from a single charge. No assembly required; just safety, comfort and convenience all at your fingertips.

KEY BENEFITS



165KG LOAD CAPACITY



TRAVEL FRIENDLY



ONLY 28KG



EASY TO USE



SMART BRAKING SYSTEM



LIGHTWEIGHT & COMPACT

Simply pull a handle and lift the seat, and watch it unfold from being compact and transportable, to a comfortable and simple-to-use solution to all your portable mobility scooter needs!

EASY TO USE

1



2



3



SPECIFICATION

BODY

- Dimensions (Unfolded):
H: 87cm W: 56cm L: 98cm
- Dimensions (Folded):
H: 54 cm W: 56cm L: 98 cm
- Ground clearance: 4.7 inches / 12cm

WHEELS

- Wheel size: diameter F: 25 cm / R: 25 cm
- Wheel type: Pneumatic
- Maximum incline: 15 degrees +-1 degree

BATTERIES

- 1 removable lithium-ion battery that goes for 30 - 35km / 18m - 21m on a single charge
- Charging time: 10 - 12 hours

KEY FEATURES

- ✓ Weighs Only 28kg
- ✓ Load Capacity: 165kg / 25 stone
- ✓ Small Enough To Fit Into The Boot Of A Car
- ✓ Safety Light

- ✓ Anti Tip Wheels
- ✓ Adjustable Handlebar
- ✓ Intelligent Braking System
- ✓ Easy Speed Control

- ✓ Removable Battery
- ✓ Expandable Footrest
- ✓ Luggage Carrier
- ✓ Easy To Fold & Store

FAQ

How do I claim back VAT?

Purchase your scooter and give us a call or drop us an email and we can refund the VAT back to you if you can provide us with a reason for your VAT relief.

How does this scooter controlled?

This scooter is operated with a button control panel found on the top of the handlebar. Ensure the scooter is on and a mode is selected to use.

Does this scooter come in other colours?

Yes this scooter does come in other colours. We offer a Red or a Black design of this model.

How do I charge the battery?

Remove the battery from the scooter and connect your charger to the port and plug the mains into the wall. Batteries will take around 10-12 hours to fully charge.

How can I put this scooter in my car?

Ensure the scooter is folded down and ensure to lift the scooter from the base as this will provide the most grip and distribute the weight evenly.

Is this scooter waterproof?

This scooter is not waterproof and we suggest drying it off with a cloth and removing the battery if in contact with heavy rain.

£2250 INC VAT

£1850 EX VAT

Check out our...

AUTO

FOLDING

MOBILITY SCOOTER

ZINNIAX



UK'S MOST CONVENIENT LIGHTWEIGHT FOLDING MOBILITY SCOOTER

Travel

Zinnia X's handlebar set contains an LCD display, so you can monitor your battery-life, speed, temperature and distance driven with ease. The lithium batteries can be removed and connected to a power outlet for charging, and offer a travel range of up to 20km.

Safety

Our intelligent braking system makes sure that whenever you're not accelerating or reversing your mobility scooter, the electromagnetic brakes are automatically applied, giving you full control wherever you decide to stop.

Convenience

This lightweight mobility scooter is incredibly travel-friendly; folding compactly for ease of transport, it weighs 32kgs and can fit snugly in the boot of a car, so you can take it anywhere! Zinnia X can go on trains, buses and coaches. And if you are looking to go a little further afield then it can also be taken on planes and cruises, making your holidays easy and accessible.



KEY BENEFITS



AIRLINE FRIENDLY



TRAVEL FRIENDLY



ONLY 32KG



EASY TO USE



SMART BRAKING SYSTEM



AUTOMATIC FOLDING

At the press of a button, watch it unfold from being compact and transportable, to a comfortable and simple-to-use solution to all your travel mobility scooter needs!

EASY TO USE

1



2



3



4



MANUAL FOLD

AUTOMATIC FOLD

SPECIFICATION

BODY

- Dimensions (Unfolded):
H: 83cm W: 49cm L: 93cm
- Dimensions (Folded):
H: 83cm W: 49cm L: 39cm
- Ground clearance: 5.2 inches / 13.5cm

WHEELS

- Wheel size: diameter F: 16cm /
R: 18cm
- Wheel type: PU wheel
- Maximum incline: : 12 degrees
+1 degree

BATTERIES

- 2 removable lithium-ion batteries
that go for 20km (12 miles) on a
single charge
- Charging time: 8 - 12 hours
- Airline friendly

KEY FEATURES

- ✓ Weighs Only 32kg
- ✓ Carrying Weight Of Up To 125kg
- ✓ Folds Small Enough To Fit Into
The Boot Of A Car
- ✓ LCD Display Unit

- ✓ Anti Tip Wheels
- ✓ Adjustable Height & Tilt
- ✓ Intelligent Braking System
- ✓ Easy Speed Control

- ✓ Removable Batteries
- ✓ Automatic Folding
- ✓ Airline Friendly
- ✓ Easy To Fold & Store

FAQ

How do I claim back VAT?

Purchase your scooter and give us a call or drop us an email and we can refund the VAT back to you if you can provide us with a reason for your VAT relief.

How does this scooter controlled?

This scooter is operated with a key fob to fold/unfold and with a dial found on the top of the LCD unit to control the speed.

Does this scooter come in other colours?

Yes this scooter does come in other colours. We offer a Red or a Blue design of this model.

How do I charge the battery?

Remove the battery from the scooter and connect your charger to the port and plug the mains into the wall. Or use the adapter to charge the batteries while in use.

How can I put this scooter in my car?

Ensure the scooter is folded down and ensure to lift the scooter from the base as this will provide the most grip and distribute the weight evenly.

Is this scooter waterproof?

This scooter is not waterproof and we suggest drying it off with a cloth and removing the battery if in contact with heavy rain.

£1950 INC VAT

£2340 EX VAT

Why not

COMPARE

OUR

PRODUCTS?

COMPARISON TABLE

COMPARE	LUPIN EASY FOLDING	REN EASY FOLDING	ZINNIA X AUTO FOLDING
WEIGHT	21KG (WITHOUT BATTERY) 23 KG (WITH BATTERY)	26 KG (WITHOUT BATTERY) 28 KG (WITHOUT BATTERY)	30KG (WITHOUT BATTERY) 32 KG (WITHOUT BATTERY)
RANGE	35KM / 21 MILES	35KM / 21 MILES	20 KM / 12.4 MILES
SPEED	4, 6, 8 MPH	4, 6, 8 MPH	0 - 3.7 MPH
LOAD CAPACITY	100 KG / 15.7 STONE	165KG / 25 STONE	125 KG / 19.6 STONE
DIMENSIONS (FOLDED)	H: 54CM W: 56CM L: 98CM	H: 54CM W: 56CM L: 98CM	H: 83CM W: 49CM L: 47CM
DIMENSIONS (UNFOLDED)	H: 87CM W: 56CM L: 98CM	H: 87CM W: 56CM L: 98CM	H: 83CM W: 49CM L: 93CM
GROUND CLEARANCE	12 CM / 4.7 INCHES	12 CM / 4.7 INCHES	13.5 CM / 5.3 INCHES
FOLDING TYPE	MANUAL	MANUAL	AUTOMATIC
BATTERY TYPE	360 Wh LITHIUM ION	374 WH LITHIUM ION	144 Wh LITHIUM ION
BATTERY QTY	1	1	2
CHARGING TIME	10 - 12 HOURS	10 - 12 HOURS	8 - 12 HOURS
ADJUSTABLE	HANDLEBAR HEIGHT	HANDLEBAR HEIGHT	ARMRESTS, HANDLEBAR HEIGHT & TILT
INCLINE	8°±1°	15°±1°	12°±1°
AIRLINE FRIENDLY	YES (WITH AIRSAFE BATTERY)	NO	YES
WHEEL SIZE (DIAMETER)	F: 25 CM / 9.8 IN B: 25 CM / 9.8 IN	F: 25 CM / 9.8 IN B: 25 CM / 9.8 IN	F: 16 CM / 6IN B: 18 CM / 7IN
DRIVE TYPE	REAR WHEEL	DUAL REAR WHEEL	REAR WHEEL
WHEEL TYPE	PNEUMATIC AIR FILLED TYRES	PNEUMATIC AIR FILLED TYRES	PU SOLID TYRES

Check out our

**ACCESSORIES
& CAR BOOT
HOISTS**

ACCESSORIES

BAGS & BASKETS



Explore our curated collection of Bags & Baskets for your mobility scooter, offering practicality and style in equal measure. From convenient storage solutions to trendy accessories, find the perfect addition to enhance your scooter experience. Stay organised and effortlessly carry your essentials, combining functionality with fashion. Discover our range now and elevate your mobility with ease.



BATTERIES & CHARGERS



Discover practical and stylish solutions for your mobility scooter with our curated collection of Bags & Baskets. From convenient storage options to trendy accessories, find the perfect addition to enhance your scooter experience. Stay organized and carry your essentials effortlessly, combining functionality with fashion. Explore our range now and elevate your mobility with ease.



COVERS & PONCHOS

Shield your mobility scooter from the elements with our premium collection of safety accessories. Explore waterproof covers and ponchos designed to protect your scooter from rain and adverse weather conditions.



SAFETY

Prioritise safety and security with our curated collection of mobility scooter accessories. From theft alarms to puncture-resistant inner tubes, our range is designed to safeguard your mobility and enhance your peace of mind.



CAR BOOT HOISTS

ATLAS 4 FOLDING MOBILITY SCOOTER HOIST



The Atlas 4 Folding Mobility Scooter allows you to easily and conveniently get your folding mobility scooter in and out of the boot of your car. Being able to get your scooter in and out of your car with ease allows you the freedom to go where you want and then enjoy it on your scooter.

MOBILITY IN MOTION

Betty and Bertie Scooters are ideal for the Autochair Lifter LM Range. With a max lifting weight of either 40kg or 80kg, this lifter is perfect for smaller, lighter mobility vehicles and can help you get out and about quickly and easily. With a fully automated lifting function, operated through the touch of a button, this lifter will eliminate the need to pick up your scooter when putting into or taking from the car.



CONTACT US

MON - FRI: 09:00 - 17:00



PHONE US
0800 999 6606



EMAIL US
info@bettyandbertie.com



VISIT US
Unit H, The Courtyard Business Centre,
Lonsome Lane, Reigate, Surrey, RH2 7QT




HOME TRIAL
Take a look online for more
information on our 14 day home trial

CHECK OUT OUR SOCIAL MEDIA TOO



FACEBOOK
[@bettyandbertiemobility](https://www.facebook.com/bettyandbertiemobility)



YOUTUBE
[@Betty&Bertie](https://www.youtube.com/Betty&Bertie)