

Product Data Sheet



Product Name *3-D Life* CD-Link

Catalog Number L60-1

Description *3-D Life* CD-Link is a component of the *3-D Life* Hydrogel system. The CD-Link molecule consists of a polyethylene glycol-peptide conjugate carrying a thiol at each end. When combined with polymers of the *3-D Life* Hydrogel system, thiol groups on CD-Link form stable thioether bonds with the thiol-reactive groups on the polymers, which results in the formation of the hydrogel. CD-Link can be used with *3-D Life* functionalised polymers Mal-PVA, Mal-Dextran, SG-PVA and SG-Dextran.

CD-Link contains a matrix metalloproteinase (MMP)-cleavable peptide sequence (Pro-Leu-Gly-Leu-Trp-Ala), which is cleaved by a broad range of MMPs including MMP1, MMP3, MMP7 and MMP9 [1]. It allows cells to locally degrade the polymer network, if they produce the indicated MMPs.

For instructions, please consult General Protocols GP-1 "Preparation of *3-D Life* Fast Gelling Hydrogels" or GP-2 "Preparation of *3-D Life* Slow Gelling Hydrogels" and the *3-D Life* Hydrogels User Guide on www.cellendes.com.

Quantity When used with polymers of the *3-D Life* Hydrogel system up to 2 ml *3-D Life* Hydrogel can be generated, depending on the stiffness of the gel.

Components

Material	Quantity	Concentration of reactive groups	Storage
● CD-Link, lyophilized	200 µl*	20 mmol/L*	Lyophilisate and after reconstitution: -20°C to -80°C
○ Water	600 µl	n/a	Room temperature or lower

All materials are filter-sterilized.

*Volume/concentration after reconstitution of lyophilisate.

Reconstitution CD-Link:

1. Briefly centrifuge vial containing CD-Link lyophilisate to make sure that the entire material is at the bottom of the reaction tube.
2. Add 188 µl *3-D Life* Water per tube for a concentration of 20 mmol/L thiol groups.
3. Close tube and briefly vortex.
4. Incubate for 5 min.
5. Briefly vortex and centrifuge again.
6. CD-Link is now ready for use.

Intended for research use only. Not for use in human therapeutic or diagnostic applications.



Distributed by: Ilex Life Sciences LLC

Tel: (828) 531-9949, Email: info@ilexlife.com, Web: <https://ilexlife.com/>