

FreeStyle® Comfort®

The perfect fit



Key features and benefits



Optimal Comfort

Let your patients experience the comfort of the body-skimming, curved silhouette design. The FreeStyle Comfort portable oxygen concentrator is lightweight and compact.



Smart O2 delivery

Unique clinical features including autoSAT®, UltraSense®, autoDOSE, and myCAIRE™ help to ensure your patients will get the oxygen they need, when they need it.



Proper O2 Saturation

With 5 settings and up to 1050 mL of oxygen per minute, the FreeStyle Comfort serves a wide range of patient needs. The unit delivers up to 95% pure oxygen to further ensure proper saturation.



Proven Quality

Robust design features, including the sieve bed technology and critical components, allow for minimum maintenance and long life. CAIRE offers one of the best warranties available, with a 3-year warranty on the device and sieve beds.

Specifications

FLOW RATE

Pulse Dose 1-5
O2 Output Settings:
Setting 1-210 mL/min
Setting 2-420 mL/min
Setting 3-630 mL/min
Setting 4-840 mL/min
Setting 5-1050 mL/min

DIMENSIONS

10.0 x 7.3 x 3.1 in
(25.4 x 18.5 x 7.9 cm)
8-cell battery pack
11.0 x 7.3 x 3.1 in
(27.9 x 18.5 x 7.9 cm)
16-cell battery pack

WEIGHT

5 lbs (2.3 kg)
8-cell battery pack
6 lbs (2.7 kg)
16-cell battery pack

SOUND LEVEL

39.9 dB(A) @ Setting 2

TRIGGER SENSITIVITY

>0.10 cm/H2O
(adjustable)

BATTERY DURATION

8-cell battery
up to 8 hours
16-cell battery
up to 16 hours

POWER

AC Power: 1-0240
VAC, 50-60Hz
DC Power: 11-18 VDC
(10 max amp)
Battery: Lithium Ion

WARRANTY

3 years
(Includes sieve beds)
1 year on battery and
power supplies



FreeStyle® Comfort®

Open the door to a bigger,
more enjoyable life

**CARRYING BAG
(Holds FreeStyle Comfort only)**

Item Number
MI396-1



**DESKTOP BATTERY CHARGER
(Battery Not Included)**

Item Numbers
BT036-1
(Type B Plug)*



BT036-2
(Type E Plug)**



**CARRY-ALL ACCESSORY BAG
(Holds Unit and Accessories)**


Item Number
MI372-2




**Replacement Strap
Item Number**
MI397-1

BATTERY PACK

8-Cell Item Number
BT034-1S



16-Cell Item Number
BT035-1S



**VEHICLE POWER CORD
(12 VDC Outlet)**

Item Number
CDO41-1




AU POWER CORD

Item Number
21347183



**POWER SUPPLY (120W 15V)
(Fits All Plug Styles)**

Item Number
PW036-1



BACKPACK

Item Number
MI459-1




TRAVEL ADAPTER

Item Number 5591-SEQ



**POWER SUPPLY WITH
AC INPUT CABLE**

Item Numbers
PW036-1S
(Type A Plug)*




PW036-2S
(Type C Plug)*



AC INPUT CABLE

Item Numbers
4997-SEQ (Type A Plug)*
4998-SEQ (Type C Plug)*



GROSS PARTICLE FILTERS

Right Filter Item Number	Left Filter Item Number
F1226-1	F1227-1



* See <https://www.iec.ch/worldplugs/map.htm#> for details on plug types.

** Purchase 5591-SEQ along with this item to properly adapt the plug for IEC type I power sockets.

Find out more about FreeStyle Comfort and accessory products at

CAIREinc.com

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.