

## Home Oxygen Concentrator



## Inogen At Home System (GS-100):

Only 8,16 kg. 5 Liters Per Minute, Continuous Flow Energy Efficient



## **Innovation in Home Oxygen Therapy**

The Inogen At Home is one of the lightest and most energy-efficient 5 liter per minute continuous flow oxygen concentrators on the market today.

At approximately 8,16 kg., it is half the weight of some other oxygen concentrators in use today. At very low power consumption levels, the Inogen At Home is an efficient choice for your next home oxygen concentrator.

Weight	8,16 kg. (18 lbs.)
Size	Width: 33,02 cm. (13.0 in.) Depth: 17,78 cm. (7.0 in.) Height: 41,91 cm. (16.5 in.)
Oxygen Flow	Continuous flow system Five flow settings from 1-5 liters per minute
Power Sensing	AC Power Supply 100-240V, 275W, 50-60Hz (auto sensing to allow worldwide use)
Operation	Simple control functions
Use	Designed for stationary use
Sound	45 dBA**

\*\*Per MDS-H1 at 2LPM

