




UVLI Slit-Top Covers 0751

Data Sheet

For 1/2-liter IV bags (500 ml). Green UVLI™ covers provide a convenient, cost-effective, and disposable solution for protecting light-sensitive medications during transportation and administration. The semi-transparency of our Green UVLI™ makes it easy to read labels and levels underneath without lifting the cover. Our covers and bags exceed the United States Pharmacopoeia Standard blocking over 90% of harmful UV light between 290nm and 450nm, so you can protect with confidence. Slit-top UVLI™ covers are fully sealed on two sides with the top seal having a convenient unsealed portion for even easier hanging.

Color	<i>Green</i>
Standard Size	<i>6 in x 10 in</i>
Legacy SKU	<i>751-SGHL</i>
Metric Size	<i>15,2 cm x 25,4 cm</i>
Quantity	<i>1000 Bags</i>
Closure	<i>Sealed on three sides (end seal)</i>
Tolerance	<i>+/- .25 in width, +/- .5 in length</i>
Capacity	<i>Not tested for weight capacity</i>
Material	<i>Low-Density Polyethylene (LDPE)</i>
FDA	<i>Approved for direct contact with polymers and colorants meeting FDA requirements and generally recognized as safe</i>
CE	<i>In Process</i>
Package Weight	<i>1 lb</i>
EPS Size	<i>Small</i>
Origin	<i>USA</i>



NOTICE: Because use conditions and applicable laws may differ from one location to another and may change with time, the Customer is responsible for determining whether products and the information in this document are appropriate for the Customer's use and for ensuring that the Customer's workplace and disposal practices are in compliance with applicable laws and other government enactments. Medipak or its affiliated company assumes no obligation or liability for the information in this document. No warranties are given; all implied warranties of merchantability or fitness for a particular purpose are expressly excluded.

Updated Date
November 16, 2022

™ Trademark of JRE Polymers, Inc. and Medipak
Supersedes all versions prior to the Updated Date