



UVLI Slit-Top Covers 0881

Data Sheet

Ideal for usage with 3-liter IV bags (3000 ml). Black UVLI™ covers provide a convenient, cost-effective, and disposable solution for protecting light-sensitive medications during transportation and administration. Our Black UVLI™ covers block 99% or more of light transmission and are opaque, making them the perfect solution for blinding. Slit-top UVLI™ covers are fully sealed on two sides with the top seal having a convenient unsealed portion for even easier hanging.

Color	<i>Black</i>
Standard Size	<i>10 in x 18 in</i>
Legacy SKU	<i>881-SB3L</i>
Metric Size	<i>25,4 cm x 45,7 cm</i>
Quantity	<i>250 Bags</i>
Closure	<i>Sealed on three sides (end seal)</i>
Tolerance	<i>+/- .25 in width, +/- .5 in length</i>
Capacity	<i>Not tested for weight capacity</i>
Material	<i>Low-Density Polyethylene (LDPE)</i>
FDA	<i>Approved for direct contact with polymers and colorants meeting FDA requirements and generally recognized as safe</i>
CE	<i>In Process</i>
Package Weight	<i>1 lb</i>
EPS Size	<i>Small</i>
Origin	<i>USA</i>



NOTICE: Because use conditions and applicable laws may differ from one location to another and may change with time, the Customer is responsible for determining whether products and the information in this document are appropriate for the Customer's use and for ensuring that the Customer's workplace and disposal practices are in compliance with applicable laws and other government enactments. Medipak or its affiliated company assumes no obligation or liability for the information in this document. No warranties are given; all implied warranties of merchantability or fitness for a particular purpose are expressly excluded.

Updated Date
November 28, 2022

™ Trademark of JRE Polymers, Inc. and Medipak
Supersedes all versions prior to the Updated Date