




UVLI Slit-Top Covers 0951

Data Sheet

For 1/2-liter IV bags (500 ml). Brown UVLI™ covers provide a convenient, cost-effective, and disposable solution for protecting light-sensitive medications during transportation and administration. Our Brown UVLI™ covers block 98% or more of light transmission and are nearly opaque, making them a great solution for blinding. Slit-top UVLI™ covers are fully sealed on two sides with the top seal having a convenient unsealed portion for even easier hanging.

Color	<i>Brown</i>
Standard Size	<i>6 in x 10 in</i>
Legacy SKU	<i>951-SDHL</i>
Metric Size	<i>15,2 cm x 25,4 cm</i>
Quantity	<i>1000 Bags</i>
Closure	<i>Sealed on three sides (end seal)</i>
Tolerance	<i>+/- .25 in width, +/- .5 in length</i>
Capacity	<i>Not tested for weight capacity</i>
Material	<i>Low-Density Polyethylene (LDPE)</i>
FDA	<i>Approved for direct contact with polymers and colorants meeting FDA requirements and generally recognized as safe</i>
CE	<i>In Process</i>
Package Weight	<i>1 lb</i>
EPS Size	<i>Small</i>
Origin	<i>USA</i>



NOTICE: Because use conditions and applicable laws may differ from one location to another and may change with time, the Customer is responsible for determining whether products and the information in this document are appropriate for the Customer's use and for ensuring that the Customer's workplace and disposal practices are in compliance with applicable laws and other government enactments. Medipak or its affiliated company assumes no obligation or liability for the information in this document. No warranties are given; all implied warranties of merchantability or fitness for a particular purpose are expressly excluded.

Updated Date
November 1, 2022

™ Trademark of JRE Polymers, Inc. and Medipak
Supersedes all versions prior to the Updated Date