




UVLI Tubing Covers 0590

Data Sheet

Convenient for use over IV Lines when cut to your desired length or for self-sealing syringes. Amber UVLI™ tubing covers provide a convenient, cost-effective, and disposable solution for protecting light-sensitive medications during administration. The semi-transparency of our Amber UVLI™ makes it easy to read labels and levels underneath without lifting the cover. Our covers and bags exceed the United States Pharmacopoeia Standard blocking over 90% of harmful UV light between 290nm and 450nm, so you can protect with confidence. Tubing covers are unsealed on either end and come in various continuous lengths for custom cutting or convenient packs of 72" precut strips for ultimate convenience.

Color	<i>Amber</i>
Standard Size	<i>2.5 in x 2500 ft</i>
Legacy SKU	<i>590-UT</i>
Metric Size	<i>6,3 cm x 762 m</i>
Quantity	<i>1 Roll</i>
Closure	<i>Sealed on three sides (end seal)</i>
Tolerance	<i>+/- .25 in width, +/- .5 in length</i>
Capacity	<i>Not tested for weight capacity</i>
Material	<i>Low-Density Polyethylene (LDPE)</i>
FDA	<i>Approved for direct contact with polymers and colorants meeting FDA requirements and generally recognized as safe</i>
CE	<i>In Process</i>
Package Weight	<i>1 lb</i>
EPS Size	<i>Small</i>
Origin	<i>USA</i>



NOTICE: Because use conditions and applicable laws may differ from one location to another and may change with time, the Customer is responsible for determining whether products and the information in this document are appropriate for the Customer's use and for ensuring that the Customer's workplace and disposal practices are in compliance with applicable laws and other government enactments. Medipak or its affiliated company assumes no obligation or liability for the information in this document. No warranties are given; all implied warranties of merchantability or fitness for a particular purpose are expressly excluded.

Updated Date
August 17, 2022

™ Trademark of JRE Polymers, Inc. and Medipak
Supersedes all versions prior to the Updated Date