




UVLI Tubing Covers 990X

Data Sheet

For IV Lines. Brown UVLI™ tubing covers provide a convenient, cost-effective, and disposable solution for protecting light-sensitive medications during administration. Our Brown UVLI™ covers block 98% or more of light transmission and are nearly opaque, making them a great solution for blinding. Tubing covers are unsealed on either end and come in various continuous lengths for custom cutting or convenient packs of 72" precut strips for ultimate convenience.

Color	<i>Brown</i>
Standard Size	<i>2.5 in x 72 in</i>
Legacy SKU	<i>990-DT6BX</i>
Metric Size	<i>6,3 cm x 1,8 m</i>
Quantity	<i>100 Precut</i>
Closure	<i>Sealed on three sides (end seal)</i>
Tolerance	<i>+/- .25 in width, +/- .5 in length</i>
Capacity	<i>Not tested for weight capacity</i>
Material	<i>Low-Density Polyethylene (LDPE)</i>
FDA	<i>Approved for direct contact with polymers and colorants meeting FDA requirements and generally recognized as safe</i>
CE	<i>In Process</i>
Package Weight	<i>1 lb</i>
EPS Size	<i>Small</i>
Origin	<i>USA</i>



NOTICE: Because use conditions and applicable laws may differ from one location to another and may change with time, the Customer is responsible for determining whether products and the information in this document are appropriate for the Customer's use and for ensuring that the Customer's workplace and disposal practices are in compliance with applicable laws and other government enactments. Medipak or its affiliated company assumes no obligation or liability for the information in this document. No warranties are given; all implied warranties of merchantability or fitness for a particular purpose are expressly excluded.

Updated Date
June 10, 2022

™ Trademark of JRE Polymers, Inc. and Medipak
Supersedes all versions prior to the Updated Date