

Turtlebrace Infant, Pediatric, and Adult Ankle Brace (AB) Sizing Chart

The following measurements are the physical dimensions before warming.

- When the brace is warmed and soft, it can be stretched more notably in circumference than length. See estimates of stretched measurements.
- The potential for stretch is proportional to the size of the brace. Larger sizes stretch further than smaller sizes
- Stretch in shaft and foot length is limited by the nonelastic zipper.
- You can expect some slight shrinkage (up to 2%) after the fitted brace cools, with less shrinking if you wait until the brace is fully hardened before removal.

INFANT	Models	Ankle Circumference		Foot Length		Shaft Height		Mid-Calf Circumference	
		Minimum + ~10 stretched		Minimum + ~3 stretched		Minimum + ~3 stretched		Minimum + ~10 stretched	
	TBBC-01 Small -1M		10.8 cm		6.5 cm		11 cm		12.6 cm
	TBBC -02 Medium - 6M		12.6 cm		8 cm		15 cm		16.2cm
	TBBC-03 Large - 18M		14.4 cm		9 cm		19 cm		19.8 cm

	Models	Ankle Circumference		Foot Lengt	h	Shaft Height	
		Minimum + ~10 stretched		Minimum + ~3 stretched		Minimum + ~3 stretched	
IATRIC	TBCP-01 Small - 5Y		16.5 cm		9 cm		22 cm
PEDI	TBCP-02 Medium - 7Y		17.8 cm		13 cm		28 cm
	TBCP-03 Large - 8Y		19 cm		13.5 cm		29 cm

	Models	Ankle Circumference		Foot Length		Shaft Height	
		Minimum + ~10 stretched		Minimum + ~3 stretched		Minimum + ~3 stretched	
ULT	TBCA-01 Small		20.3cm		18 cm		32 cm
AD	TBCA-02 Medium		21.6cm		19 cm		36 cm
	TBCA-03 Large		22.9cm		21.5 cm		39 cm