Compatibility With Eye Drops



Manufacturer	Brand	Bottle Size	Opticare	Arthro 5	Arthro 10
Allergan	Combigan	5ml	\Diamond		\Diamond
	Alphagan	5ml			\Diamond
	Lumigan	3ml	\Diamond	\Diamond	
	Betagan	5ml			
	Propine	5&10ml	\Diamond		
	Ganfort	5ml			
	Pred Forte	5ml			
	Pred Forte	10ml			
	FML	5ml	\Diamond		
	FML	10ml			\Diamond
	Exocin	5ml	\Diamond		
	Acular	5ml			
	Liquifilm Tears	15ml	\Diamond		\Diamond
	Refresh Contacts	15ml			

Manufacturer	Brand	Bottle Size	Opticare	Arthro 5	Arthro 10
Bausch & Lomb	Artelac	10ml	\Diamond		
	Floxal	5ml		\Diamond	
	Betamann 0,1%	5ml		\Diamond	
	Betamann 0,3%	5ml	\Diamond		
	Timomann 0,25%	5ml	\Diamond		
Martindale	Hypromellose	10ml	\Diamond		
Pharmaceuticals	Timolol	5ml			
	Chloramphenicol	10ml			
	Homatropine	10ml			
	Neomycin Sulphate	10ml			
	Phenylephrine	10ml	\Diamond		
	Sodium Chloride	10ml			
	Vista Methasone	5ml	\Diamond	\Diamond	
	Vista Methasone	10ml			\Diamond
	Zinc Sulphate	10ml			\Diamond

Manufacturer	Brand	Bottle Size	Opticare	Arthro 5	Arthro 10
Alcon	Azopt	5ml	\Diamond	\Diamond	
	Betoptic 0,5%	5ml			
	Betoptic 0,25%	5ml	\Diamond		
	lopidine	5ml			
	Travatan #	2.5ml	\Diamond		
	Duotrav #	2.5ml			
MSD	Timoptol *	5ml	\Diamond		
	Timoptol LA **	5ml			
	Trusopt *	5ml	\Diamond		
	Cosopt *	5ml			
Pfizer	Xalatan #	2.5ml			
	Xalacom #	2.5ml			

Use bottle collar adaptor (supplied with Opticare) to ensure secure fit inside Opticare.

^{*} This bottle is more rigid than many. Insert bottle into Opticare with the top removed and flat side first.

^{**} Aid available from MSD.

Opticare eye drop dispensers are compatible with traditionally shaped bottles.

Our aids cannot be used with multi-dose preservative-free bottles with a wide neck.



