












Symbol Glossary

SYMBOL	SYMBOL TITLE	SYMBOL DESCRIPTION	SOURCE / REFERENCE
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1 #5.1.5 ISO 7000-2492
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1 #5.1.6 ISO 7000-2493
	Date of manufacture	Indicates the date when the medical device was manufactured.	ISO 15223-1 #5.1.3 ISO 7000-2497
	Manufacturer	Indicates the medical device manufacturer.	ISO 15223-1 #5.1.1 ISO 7000-3082
	Serial number	Indicates the manufacturer's serial number so that a specific medical device can be identified.	ISO 15223-1 #5.1.7 ISO 7000-2498
	Unique Device Identifier	Indicates a carrier that contains unique device identifier information.	ISO 15223-1 #5.7.10
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1 #5.4.3 ISO 7000-1641
implantwisedental.com/symbolglossary implantwisedental.com/ifu	Link to Online Symbols Glossary and IFU Portal	Identifies the web address where electronic Instructions for Use (IFU) and symbols glossary can be obtained.	ImplantWise custom symbol
	Patient information website	Indicates a website where a patient may obtain additional information on the medical product.	ISO 15223-1 #5.7.4 ISO 7000-3705
	Caution	Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.	ISO 15223-1 #5.4.4 ISO 7000-0434A
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1 #5.4.2 ISO 7000-1051
	Do not use if package is damaged and consult instructions for use	Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information.	ISO 15223-1 #5.2.8 ISO 7000-2606
	Non-sterile	Indicates a medical device that has not been subjected to a sterilization process.	ISO 15223-1 #5.2.7 ISO 7000-2609
	Medical device	Indicates the item is a medical device.	ISO 15223-1 #5.7.7

SYMBOL	SYMBOL TITLE	SYMBOL DESCRIPTION	SOURCE / REFERENCE
Rx only	For prescription use only	Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed dental practitioner.	FDA Final rule, Use of symbols in labeling, 15 June 2016