


[HTTPS://WITHOMED.COM/](https://withomed.com/)

iE 10

Mobile ECG



iE 10

Mobile ECG



Accuracy and reliability in ECG interpretation

- Biocare CardioPro™ ECG analysis program, proven by CSE/AHA/MIT database
- 12-lead simultaneous acquisition
- Complete digital filtering and self-adapting filtering
- Up to 300 seconds of rhythm analysis for easier in arrhythmia locating
- Applied the latest Minnesota code

Streamlined workflow, new operating experience

- Onscreen operation guide to check the ECG procedure
- Multi-language voice reminder and visual alarming for lead-off, and low battery
- Support general settings according to different situations for physicians and patients
- 3 report formats: simple report, median beat report and detailed report
- Support external laser printer

Data management solutions

- Massive local storage capacity up to 3000 files
- Upload the ECG data to the ECG management software via WiFi or 4G network



- About 0.5kg, light in weight



- High brightness 7.0-inch capacitive touch screen, 1280*800



- Horizontal or vertical screen selectable



- Built-in rechargeable battery, up to 4.0 hours of continuous work



WITHOMED

Add: Longhua District, Shenzhen 1301, Building 3, No. 1301-80, Sightseeing Road, Xinlan Community, Guanlan Street Shenzhen, 518110 Guangdong, China

E-mail: contact@withomed.com

Website: <https://withomed.com/>

