



Search all content

The Exchange > Sensei Cloud > SENSEI CLOUD > Phoenix Weekly Release 0.0.0.1193



MOD kevin_moloney

Moderator



10-30-2017 01:57 PM

Phoenix Weekly Release 0.0.0.1193

Phoenix Online DPMS

Weekly Release

Build 0.0.0.1193 | October 30, 2017

These are the highlights of the new features available in this week's Phoenix DPMS release.

Restorative Charting

This week's release was focused on preparing the hard tissue chart for future workflow support as well as ironing out some of the kinks in the current mechanics of anatomy selection. For example, we implemented several improvements to the logic for selecting anatomy within a single tooth to optimize the charting mechanism and avoid duplicate or illogical charting scenarios (e.g., inclusive anatomical selections). So, selection of specific crown surfaces on a tooth will replace selection of the entire crown (and vice versa). Similarly, selection of specific tooth anatomy will replace selection of the entire tooth (and vice versa).

Carestream DENTAL Checkout Queue (6): S. Hardin 08:00 am A. Sant 09:00 am

Smiley Dental Surgery + Terrence Molarnator x Kevin M

Home Patient Record Insurance Clinical Financials Trea

Kevin Moloney (404) 770-6780 (Self) 37 yrs (0) +

Current Conditions
No Current Conditions
Last Charted:

Charting

Current Selection Clear

Tooth: 1.2
6 surfaces selected

Chart Findings

Decay
Damaged

History

Additionally, improvements were made selection of anatomy across teeth so that the suggested / available findings to chart are more contextually-appropriate. This includes the ability to select variable tooth anatomy across the arches. These enhancements will help users avoid getting themselves into a state in which there is nothing logical to chart against the current set of selected anatomy and will become important as we implement support for charting of restorations and appliances.

Some significant refactoring of the hard tissue chart implementation was also completed this week. While much of this is invisible to the user, it is necessary for the support of future workflows such as charting restorations, appliances, and treatments while also visually representing all of these charting history events on the hard tissue chart. One more immediate effect of this work is that the current anatomy selections are much more visible, even when charting history already exists on the anatomy.

Additionally, we have improved the granularity of the age-based default charts, so that the process of setting up initial charts for young patients will likely require less clicks. Now, the default hard tissue chart created for a new patient record more accurately reflects the age-based development models so that dentition profile present in a 6.5 year old will be different than a 9.5 year old and so on. This should be of particular interest to pedodontists and to general dentists that see a wide age range of patients, which need workflow support for efficient setup of initial hard tissue charts.

The screenshot displays the Carestream Dental software interface. At the top, there is a navigation bar with tabs for Home, Patient Record, Insurance, Clinical (selected), Financials, Treatment Plan, and Files. Below this, a patient information section shows the name Michaela Gent, gender female, and contact information. The main area is divided into sections for Charting, Current Selection, Chart Findings, and History. The Charting section features a grid of tooth models and diagrams, labeled with numbers and letters: 1, 2, 3, A, B, C, D, E, 32, 31, 30, T, S, R, Q, and 25. Each tooth is represented by a 3D model and a corresponding 2D diagram.

Carestream DENTAL Checkout Queue (6): S. Hardin 08:00 am A. Santiago 09:00 am V. Simon 09:45 am J. Hutchinson 11:15 am D. Bradl 08:00 am

Smiley Dental Surgery + Terrence Molarnator x Kevin Moloney x Nye Moloney x Michaela Ger

Home Patient Record Insurance Clinical Financials Treatment Plan Files

Mickey Gent (404) 555-8787 9y7m Michael Gent (Parent) (0) +

Current Conditions: No Current Conditions Last Charted:

Current Restorations: No Current Restorations

Charting

Current Selection Clear

Chart Findings

History

Finally, some enhancements were also made to the charting history (notes) that document the clinical findings that users chart during an exam. These updates will become important as we start to implement charting of completed procedures and existing restorations and appliances. We will begin to support charting for existing restorations and completed work soon.

Claim Processing Enhancements

This week's release also focused on continuing the integration with the eClaims service and the Carestream Dental clearinghouse. As mentioned in last week's Release Notes, the claim file is currently being generated and submitted to a test account at the Carestream Dental clearinghouse and has passed initial testing. We are continuing to test submission of claim batches from Phoenix to the clearinghouse for processing. Additionally, some code refactoring was done to address an issue with refreshing the Claims List UIs. Next steps will include the implementation of insurance estimation and the handling of claim status reports back from the Carestream Dental clearinghouse.

Optimizations and Fixes

The following fixes and optimizations have been also been added to Phoenix this week.

1. Fix to the Patient Homepage to display notes in chronological order.
2. Addressed an issue to ensure that updates to the Patient Record are automatically reflected in open Patient Tabs.
3. Optimized the Patient Homepage to display UI components cleanly at 1366 x 768 resolution.
4. Fixed some bugs related to canceling the creation of a new patient record.
5. Optimized the display of contextual information in the Appointment Card, so that the patient's treatment plan items are automatically displayed when linking from a treatment plan booking context.
6. Addressed an issue with saving new notes on the Treatment Plan List UI.
7. Fixed some bugs that were interfering with the create new patient record workflows.

 Add tags



Reply

Sort: Newest to Oldest

[All forum topics](#) [< Previous Topic](#) [Next Topic >](#)



Reply...

Post Reply

Powered by
Khoros 

CORPORATE HEADQUARTERS

Carestream Dental LLC
3625 Cumberland Blvd. Ste. 700
Atlanta, GA 30339

CONTACT

[Contact Us](#)

[Support Locator](#)

[Training](#)

ABOUT US

[Blog](#)

[Media Room](#)

QUICK LINKS

[Carestream Dental Institute](#)

[Carestream Dental Website](#)



[Global Data Privacy Shield](#)

[Terms and Conditions](#)

[Privacy Policy](#)

© 2019 Carestream Dental, LLC. All Rights Reserved