



## TECHNICAL DATA SHEET Disinfectant Wipes

### Intended Use

Used to clean or disinfect hard, non-porous surface in dental offices, hospitals and other healthcare settings.

### Active Ingredients

Ethyl Alcohol 72.50 % w/w

Didecyl Dimethyl Ammonium Chloride 0.33 % w/w

Other Ingredients 27.17 % w/w

Characteristic	Specification		
Description	Soft non-woven wet wipe in perforated roll style (Refer to Appendix 1 for Specific Attribute)		
Material composition	Wipe	Polypropylene	
	Canister	HDPE	
	Cap	Polypropylene	
	Seal	Film	
	Label	Film	
Dimension	Size	Width (mm)	Length (mm)
	6" x 6.75"	152	172
	9" x 12"	229	305

### Visual



### Additional Information

Storage Condition	Shelf Life	Country of Origin
Protect from extreme cold or heat	2 years	USA

### Appendix 1 List of Item Codes and Attributes

## Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Treatment
40060C	Medicom Pro-Surface+ Disinfectant Wipe	160 Wipes/Canister, 12 Canisters/Case, 1920 Wipes/Case	6" x 6.75"	White wipes	Non-Sterile
40061C	Medicom Pro-Surface+ Disinfectant Wipe	65 Wipes/Canister, 12 Canisters/Case, 780 Wipes/Case	9" x 12"	White wipes	Non-Sterile

## Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Treatment
40060	Medicom Pro-Surface Disinfectant Wipe	160 Wipes/Canister, 12 Canisters/Case, 1920 Wipes/Case	6" x 6.75"	White wipes	Non-Sterile
40061	Medicom Pro-Surface Disinfectant Wipe	65 Wipes/Canister, 12 Canisters/Case, 780 Wipes/Case	9" x 12"	White wipes	Non-Sterile