



## TECHNICAL DATA SHEET

### Non-Woven Sponge

#### Intended Use

A hydrophilic non-sterile sponge intended to cover a wound, such as to be placed directly on a patient's wound to absorb exudates.

Characteristic	Specification		
Material Composition	50% Viscose and 50% Polyester (Latex free)		
Tensile Strength	MD	≥ 55 N/50 mm	
	CD	≥ 25 N/50 mm	
Description	Soft non-woven material, folded.		
Ply	4 ply		
Color	White		
Dimension	Size	Width	Length
	2" x 2" (5 x 5 cm)	2" (5,1 cm)	2" (5,1 cm)
	4" x 4" (10 x 10 cm)	4" (10,2 cm)	4" (10,2 cm)

#### Additional Information

Storage Condition	Country of Origin
Store under cool, dry conditions	United States

#### Appendix 1 List of Item Codes and Attributes

## Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Treatment	Marketed in USA as	Marketed in Canada as
4522-5	Medicom SafeGauze Non-Woven Sponges	200 Sponges/Sleeve, 25 Sleeves/Case, 5000 Sponges/Case Sleeve made of Kraft Paper, Shipper case is made of corrugated fiberboard.	2" x 2" (5 x 5 cm)	Non-Sterile	Medical Device Class 1	Medical Device Class 1
4544	Medicom SafeGauze Non-Woven Sponges	200 Sponges/Sleeve, 10 Sleeves/Case, 2000 Sponges/Case Sleeve made of Kraft Paper, Shipper case is made of corrugated fiberboard.	4" x 4" (10 x 10 cm)	Non-Sterile	Medical Device Class 1	Medical Device Class 1
4522	Medicom SafeGauze Non-Woven Sponges	200 Sponges/Sleeve, 20 Sleeves/Case, 4000 Sponges/Case Sleeve made of Kraft Paper, Shipper case is made of corrugated fiberboard.	2" x 2" (5 x 5 cm)	Non-Sterile	Medical Device Class 1	Medical Device Class 1