

TECHNICAL DATA SHEET


Shoe Cover

Intended Use

Non-Sterile Shoe Covers intended to be worn by operating room personnel during various procedures to protect patient/materials and the operating personnel from transfer of microorganisms, particulates, etc. Single use.

Characteristic		Specification	
Material Composition		Spunbonded Polypropylene, Latex free	
Description		Rectangular piece of non-woven fabric folded on the length and closed on each side. Features: An elastic closure on the top opening and an elastic fixed inside the sole, on the fold, covered by a piece fabric. (Refer to Appendix 1 specific attribute)	
Part of Shoe Cover		Dimension (cm)	
		Regular	X-Large
A	Shoe Cover Stretched Length	42	42
B	Shoe Cover Height	18	19
C	Opening Natural Length	12,5	12,5
C'	Opening Stretched Length	38	38
D	Ultrasonic Print	2,5	2,5

Measurement Diagram



Product Image



Additional Information

Storage Condition	Country of Origin
Store in a cool dry place	Canada

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
8006	Medicom SafeBasics Non-Conductive Shoe Covers	100 Shoe Covers/Box, 3 Boxes/Case, 300 Shoe Covers/Case	X-Large	Blue	Non-Skid	Non-Sterile
8007	Medicom SafeBasics Non-Conductive Shoe Covers	100 Shoe Covers/Box, 3 Boxes/Case, 300 Shoe Covers/Case	Regular	Blue	Non-Skid	Non-Sterile
8008	Medicom SafeBasics Non-Conductive Shoe Covers	100 Shoe Covers/Box, 3 Boxes/Case, 300 Shoe Covers/Case	Regular	Blue	Regular Sole	Non-Sterile
8008-XL	Medicom SafeBasics Non-Conductive Shoe Covers	100 Shoe Covers/Box, 3 Boxes/Case, 300 Shoe Covers/Case	X-Large	Blue	Regular Sole	Non-Sterile
8036	Medicom SafeBasics Non-Conductive Shoe Covers	100 Shoe Covers/Bag, 10 Bags/Case, 1000 Shoe Covers/Case	X-Large	White	Non-Skid	Non-Sterile