

## TECHNICAL DATA SHEET


### Bib

#### Intended Use

Bibs are made of paper intended to absorb saliva or moisture from the oral cavity during procedures. Latex-free.

Characteristic	Specification		
Material Composition	Tissue Paper & Polyethylene Film.		
Description	Folded piece of paper with a plastic backing. Embossed horizontal pattern on the body of the bib and vertical embossed pattern on the vertical borders. (Refer to Appendix 1 for specific attribute)		
Dimension	Unfolded Size	Unfolded Width (cm)	Unfolded Length (cm)
	13" x 18" (33 x 45 cm)	33,0	45,7

#### Measurement Diagram



#### Visual



#### Additional Information

Storage Condition	Country of Origin
Store in a cool dry place	China

#### Appendix 1 List of Item Codes and Attributes



## Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Treatment
8182	Medicom SafeBasics Dry-Back Bibs	800 Bibs/Case	13" x 18" (33 x 45 cm)	Blue	Non-Sterile
8183	Medicom SafeBasics Dry-Back Bibs	800 Bibs/Case	13" x 18" (33 x 45 cm)	White	Non-Sterile
8185	Medicom SafeBasics Dry-Back Bibs	800 Bibs/Case	13" x 18" (33 x 45 cm)	Green	Non-Sterile
8186	Medicom SafeBasics Dry-Back Bibs	800 Bibs/Case	13" x 18" (33 x 45 cm)	Pink	Non-Sterile