


TECHNICAL DATA SHEET
Bib: 2 Ply Tissue + 1 Ply Poly

Intended Use

Bibs are made of paper intended to absorb saliva or moisture from the oral cavity during procedures. Latex-free.

Characteristic	Specification		
Material Composition	Tissue Paper & Polyethylene Film		
Description	Folded piece of paper with a plastic backing. Embossed pattern on the body and the borders.		
Color	Refer to Appendix 1		
Folded Dimension	Folded Size (cm)	Folded Width (cm)	Folded Length (cm)
	11,5 x 16,5	11,5	16,5
Unfolded Dimension	Unfolded Size (cm)	Unfolded Width (cm)	Unfolded Length (cm)
	33 x 45	33	45
Piece Weight	2 Ply Tissue + 1 Ply Poly	7.2 g	



Additional Information

Storage Condition	Country of Origin
Store in a cool dry place	China

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Treatment
HE-PBB	Hedy Elements Patient Bibs	125 Bibs/Bag, 4 Bags/Case, 500 Bibs/Case	13"x18"	Blue	Non-Sterile
HE-PBG	Hedy Elements Patient Bibs	125 Bibs/Bag, 4 Bags/Case, 500 Bibs/Case	13"x18"	Green	Non-Sterile
HE-PBL	Hedy Elements Patient Bibs	125 Bibs/Bag, 4 Bags/Case, 500 Bibs/Case	13"x18"	Lavender	Non-Sterile
HE-PBM	Hedy Elements Patient Bibs	125 Bibs/Bag, 4 Bags/Case, 500 Bibs/Case	13"x18"	Mauve	Non-Sterile
HE-PBW	Hedy Elements Patient Bibs	125 Bibs/Bag, 4 Bags/Case, 500 Bibs/Case	13"x18"	White	Non-Sterile