



core3d
centres®



3Shape D900/900L®

Laboratory Scanning

Core3dcentres® are an authorised distributor and trainer for the 3Shape D900/900L® - the world's leading dental scanner and software supplier. Four high-res 5MP cameras and new blue LED technology provide extraordinary speed and highly accurate colour scanning. Perfect for large high volume, productive orientated labs.

Indications:

- Anterior and posterior single crowns
- Inlays
- Onlays
- Veneers
- Bridges
- Abutments
- Implant Bars

Advantages:

- 4 cameras, 5 mega pixels
- Light technology - blue LED
- Accuracy Crown & Bridge - 7 microns
- Accuracy Implant Bars - 8 microns
- Gypsum model scanning
- Impression scanning (Add-on)
- Expanded Interior Space (D900L model)
- Multi-Die scanning
- Coloured texture scanning
- Single die scan time - 15 sec - 19 sec
- 3 Unit bridge scan time - 50 sec - 70 sec



3Shape D900/900L[®]

Laboratory Scanning

System configuration

3D Scanner

D900/900L 3D Scanner specialised for dental applications

Desktop PC

High performance PC (model subject to change)

Operating system and CAD software package pre-installed

Mouse and keyboard included (Monitor NOT included)

Operating System

Windows 8, 64-bit, English (Support for 35 languages)

CAD Software package

Dental System[™] Premium, consists of:

Dental Designer

(CAD software for dental restoration, 1 client point)

Dental Manager

(Order Management software, 5 client points)

Dental Manager Inbox

(Order reception and filtering software, 3 client points)

Teamviewer Client

(Remote assistance software)

Abutment Designer[™] / Implant Bars and Bridge Design (Limited saves: 20 + 5)

FREE trial software modules

Removable Partial Design (Limited saves: 15)

Help System Integrated online user manual (11 languages)

Repair Disk

Windows 8 repair disk

Calibration Box

1 x Scanner calibration object

3 x Scanning interface plate

1 x Blu-Tack[®] scan clay

1 x Dongle USB key

Annual Licence Fee

Optional

Upgrade to laptop

High-performance laptop (model subject to change) replacing standard desktop PC. Includes external mouse

Additional Client Points

Additional client points can be purchased

Add-on Modules

Various add-on modules are available. Visit www.core3dcentres.com/solutions/software to find out more

For the laboratory

What to send:

To receive a restoration from Core3dcentres, the laboratory will be required to send a digitally designed file, which we will configure your system to allow you to send directly to our 3Shape system, with reference to the shade required and any special requirements and clear instructions of what has to be manufactured.

For more detailed information please contact us at www.core3dcentres.com.

Asia Pacific Benelux Canada Germany Japan Scandinavia South Europe Switzerland United Kingdom USA