



# 3Shape Dental System™

## Product Catalog

# 3Shape D500

D500 3D scanner with Dental System™ Full Contour



Item Number

**20001700**

20001705 - Demo system (see note 1)

Send your order to [order@3shape.com](mailto:order@3shape.com)

## System Configuration

### 3D Scanner

D500 3D Scanner (2 x 1.3 MP cameras)  
optimized for dental applications.

### Desktop PC

- High performance PC (model subject to change)
- Operating system and CAD software pre-installed
- Mouse and keyboard included (monitor NOT included)

### Operating System

- Windows 8 Pro, 64-bit, English  
(Support for 35 languages)

### Repair Disk

- Windows 8 Pro repair disk

### Accessories

- 1 x Scanner calibration object
- 3 x Scanning interface plate
- 1 x Blu-Tack® scan clay
- 1 x Dongle USB key

### CAD Software Package

Dental System™ Full Contour consisting of:

- Dental Designer - Full Contour; see page 19  
(CAD software for dental restorations, 1 seat)
- Dental Manager  
(Order Management software, 5 seats)
- Dental Manager Inbox  
(Order reception and filtering software, 5 connections)
- TRIOS® Inbox  
(TRIOS® order reception and filtering software)

### FREE Trial Software Modules

- Abutment Designer  
(Limited saves: 20)
- Removable Partial Design  
(Limited saves: 15)
- Model Builder™ - TRIOS®  
(Limited saves: 20)

### Help System

- Integrated online user manual  
(12 languages)

### Optional:

#### Add-on Modules and Module Upgrades

#### CAD Points

#### Additional Seats

#### Extended Warranty

Various add-on modules and Module Upgrades are available. See page 11-13.  
Pay-Per-Design option for selected add-on modules. Purchase in webshop.  
Additional seats can be purchased. See separate list.  
Extend the 1-year limited warranty on 3Shape hardware by 1 or 2 additional years.

#### Notes:

1. Cannot be resold for a period of 12 months. Limited to 200 saves yearly.

# 3Shape D500

D500 3D scanner with Dental System™ Standard



Item Number  
**20001600**

20001605 - Demo system (see note 2)

Send your order to [order@3shape.com](mailto:order@3shape.com)

## System Configuration

### 3D Scanner

D500 3D Scanner (2 x 1.3 MP cameras)  
optimized for dental applications.

### Desktop PC

- High performance PC (model subject to change)
- Operating system and CAD software pre-installed
- Mouse and keyboard included (monitor NOT included)

### Operating System

- Windows 8 Pro, 64-bit, English  
(Support for 35 languages)

### Repair Disk

- Windows 8 Pro repair disk

### Accessories

- 1 x Scanner calibration object
- 3 x Scanning interface plate
- 1 x Blu-Tack® scan clay
- 1 x Dongle USB key

### CAD Software Package

Dental System™ Standard consisting of:

- Dental Designer - Standard; see page 19  
(CAD software for dental restorations, 1 seat)
- Dental Manager  
(Order Management software, 5 seats)
- Dental Manager Inbox  
(Order reception and filtering software, 5 connections)
- TRIOS® Inbox  
(TRIOS® order reception and filtering software)

### FREE Trial Software Modules

- Abutment Designer  
(Limited saves: 20)
- Removable Partial Design  
(Limited saves: 15)
- Model Builder™ - TRIOS®  
(Limited saves: 20)

### Help System

- Integrated online user manual  
(12 languages)

### Optional:

**Add-on Modules and Module Upgrades**  
**CAD Points**  
**Additional Seats**  
**Extended Warranty**

Various add-on modules and Module Upgrades are available. See page 11-13.  
Pay-Per-Design option for selected add-on modules. Purchase in webshop.  
Additional seats can be purchased. See separate list.  
Extend the 1-year limited warranty on 3Shape hardware by 1 or 2 additional years.

#### Notes:

2. Cannot be resold for a period of 12 months. Limited to 200 saves yearly.

# 3Shape D700

D700 3D scanner with Dental System™ Premium



Item Number  
**20001000**

20001005 - Demo system (see note 3)  
20001010 - Educational system (see note 4)

Send your order to [order@3shape.com](mailto:order@3shape.com)

## System Configuration

### 3D Scanner

D700 3D Scanner (2 x 1.3 MP cameras)  
optimized for dental applications.

### Desktop PC

- High performance PC (model subject to change)
- Operating system and CAD software pre-installed
- Mouse and keyboard included (monitor NOT included)

### Operating System

- Windows 8 Pro, 64-bit, English  
(Support for 35 languages)

### Repair Disk

- Windows 8 Pro repair disk

### Accessories

- 1 x Scanner calibration object
- 3 x Scanning interface plate
- 1 x Blu-Tack® scan clay
- 1 x Dongle USB key

### CAD Software Package

Dental System™ Premium consisting of:

- Dental Designer - Premium; see page 19  
(CAD software for dental restorations, 1 seat)
- Dental Manager  
(Order Management software, 5 seats)
- Dental Manager Inbox  
(Order reception and filtering software, 10 connections)
- TRIOS® Inbox  
(TRIOS® order reception and filtering software)

### FREE Trial Software Modules

- Abutment Designer / Implant Bars and Bridge Design  
(Limited saves: 20/5)
- Removable Partial Design  
(Limited saves: 15)
- Model Builder™ - TRIOS®  
(Limited saves: 20)

### Help System

- Integrated online user manual  
(12 languages)

### Optional:

**Add-on Modules and Module Upgrades**  
**CAD Points**  
**Additional Seats**  
**Extended Warranty**

Various add-on modules and Module Upgrades are available. See page 11-13.  
Pay-Per-Design option for selected add-on modules. Purchase in webshop.  
Additional seats can be purchased. See separate list.  
Extend the 1-year limited warranty on 3Shape hardware by 1 or 2 additional years.

#### Notes:

3. Cannot be resold for a period of 12 months. Limited to 200 saves yearly.

4. Licenses are renewed upon request every year. All upgrades are included for a period of 3 years. The name of the school or university and name of the contact person are required at order placement.

# 3Shape D710

D710 3D scanner with Dental System™ Premium



Item Number

**20001100**

20001105 - Demo system (see note 5)

20001110 - Educational system (see note 6)

Send your order to [order@3shape.com](mailto:order@3shape.com)

## System Configuration

### 3D Scanner

D710 3D Scanner (2 x 1.3 MP cameras)  
optimized for dental applications. With Multi-Die Scanning.

### Desktop PC

- High performance PC (model subject to change)
- Operating system and CAD software pre-installed
- Mouse and keyboard included (monitor NOT included)

### Operating System

- Windows 8 Pro, 64-bit, English  
(Support for 35 languages)

### Repair Disk

- Windows 8 Pro repair disk

### Accessories

- 1 x Scanner calibration object
- 3 x Scanning interface plate
- 1 x Blu-Tack® scan clay
- 1 x Dongle USB key
- 1 x D710/810 Multi-Die Fixture

### CAD Software Package

Dental System™ Premium consisting of:

- Dental Designer - Premium; see page 19  
(CAD software for dental restorations, 1 seat)
- Dental Manager  
(Order Management software, 5 seats)
- Dental Manager Inbox  
(Order reception and filtering software, 10 connections)
- TRIOS® Inbox  
(TRIOS® order reception and filtering software)

### FREE Trial Software Modules

- Abutment Designer / Implant Bars and Bridge Design  
(Limited saves: 20/5)
- Removable Partial Design  
(Limited saves: 15)
- Model Builder™ - TRIOS®  
(Limited saves: 20)

### Help System

- Integrated online user manual  
(12 languages)

### Optional:

#### Add-on Modules and Module Upgrades

#### CAD Points

#### Additional Seats

#### Extended Warranty

Various add-on modules and Module Upgrades are available. See page 11-13.  
Pay-Per-Design option for selected add-on modules. Purchase in webshop.  
Additional seats can be purchased. See separate list.  
Extend the 1-year limited warranty on 3Shape hardware by 1 or 2 additional years.

#### Notes:

5. Cannot be resold for a period of 12 months. Limited to 200 saves yearly.

6. Licenses are renewed upon request every year. All upgrades are included for a period of 3 years. The name of the school or university and name of the contact person are required at order placement.



# 3Shape D800

D800 3D scanner with Dental System™ Premium



Item Number

**20001200**

20001205 - Demo system (see note 7)

20001210 - Educational system (see note 8)

Send your order to [order@3shape.com](mailto:order@3shape.com)

## System Configuration

### 3D Scanner

D800 3D Scanner (2 x 5 MP cameras)  
optimized for dental applications. With Texture Scanning.

### Desktop PC

- High performance PC (model subject to change)
- Operating system and CAD software pre-installed
- Mouse and keyboard included (monitor NOT included)

### Operating System

- Windows 8 Pro, 64-bit, English  
(Support for 35 languages)

### Repair Disk

- Windows 8 Pro repair disk

### Accessories

- 1 x Scanner calibration object
- 3 x Scanning interface plate
- 1 x Blu-Tack® scan clay
- 1 x Dongle USB key

### CAD Software Package

Dental System™ Premium consisting of:

- Dental Designer - Premium; see page 19  
(CAD software for dental restorations, 1 seat)
- Dental Manager  
(Order Management software, 5 seats)
- Dental Manager Inbox  
(Order reception and filtering software, 10 connections)
- TRIOS® Inbox  
(TRIOS® order reception and filtering software)

### FREE Trial Software Modules

- Abutment Designer / Implant Bars and Bridge Design  
(Limited saves: 20/5)
- Removable Partial Design  
(Limited saves: 15)
- Model Builder™ - TRIOS®  
(Limited saves: 20)

### Help System

- Integrated online user manual  
(12 languages)

### Optional:

#### Add-on Modules and Module Upgrades

#### CAD Points

#### Additional Seats

#### Extended Warranty

Various add-on modules and Module Upgrades are available. See page 11-13.  
Pay-Per-Design option for selected add-on modules. Purchase in webshop.  
Additional seats can be purchased. See separate list.  
Extend the 1-year limited warranty on 3Shape hardware by 1 or 2 additional years.


### Notes:

7. Cannot be resold for a period of 12 months. Limited to 200 saves yearly.

8. Licenses are renewed upon request every year. All upgrades are included for a period of 3 years. The name of the school or university and name of the contact person are required at order placement.

# 3Shape D810

D810 3D scanner with Dental System™ Premium



Item Number  
**20001300**

20001305 - Demo system (see note 9)

20001310 - Educational system (see note 10)

Send your order to [order@3shape.com](mailto:order@3shape.com)

## System Configuration

### 3D Scanner

D810 3D Scanner (2 x 5 MP cameras)  
optimized for dental applications.

With Multi-Die Scanning and Texture Scanning.

### Desktop PC

- High performance PC (model subject to change)
- Operating system and CAD software pre-installed
- Mouse and keyboard included (monitor NOT included)

### Operating System

- Windows 8 Pro, 64-bit, English  
(Support for 35 languages)

### Repair Disk

- Windows 8 Pro repair disk

### Accessories

- 1 x Scanner calibration object
- 3 x Scanning interface plate
- 1 x Blu-Tack® scan clay
- 1 x Dongle USB key
- 1 x D710/810 Multi-Die Fixture

### CAD Software Package

Dental System™ Premium consisting of:

- Dental Designer - Premium; see page 19  
(CAD software for dental restorations, 1 seat)
- Dental Manager  
(Order Management software, 5 seats)
- Dental Manager Inbox  
(Order reception and filtering software, 10 connections)
- TRIOS® Inbox  
(TRIOS® order reception and filtering software)

### FREE Trial Software Modules

- Abutment Designer / Implant Bars and Bridge Design  
(Limited saves: 20/5)
- Removable Partial Design  
(Limited saves: 15)
- Model Builder™ - TRIOS®  
(Limited saves: 20)

### Help System

- Integrated online user manual  
(12 languages)

## Optional:

### Add-on Modules and Module Upgrades

#### CAD Points

#### Additional Seats

#### Extended Warranty

Various add-on modules and Module Upgrades are available. See page 11-13.

Pay-Per-Design option for selected add-on modules. Purchase in webshop.

Additional seats can be purchased. See separate list.

Extend the 1-year limited warranty on 3Shape hardware by 1 or 2 additional years.

### Notes:

9. Cannot be resold for a period of 12 months. Limited to 200 saves yearly.

10. Licenses are renewed upon request every year. All upgrades are included for a period of 3 years. The name of the school or university and name of the contact person are required at order placement.

# 3Shape D900

D900 3D scanner with Dental System™ Premium



Item Number

**20001400**

20001405 - Demo system (see note 11)

20001410 - Educational system (see note 12)

Send your order to [order@3shape.com](mailto:order@3shape.com)

## System Configuration

### 3D Scanner

D900 3D Scanner (4 x 5 MP cameras)  
optimized for dental applications.

With Multi-Die Scanning and Color Texture Scanning.

### Desktop PC

- High performance PC (model subject to change)
- Operating system and CAD software pre-installed
- Mouse and keyboard included (monitor NOT included)

### Operating System

- Windows 8 Pro, 64-bit, English  
(Support for 35 languages)

### Repair Disk

- Windows 8 Pro repair disk

### Accessories

- 1 x Scanner calibration object
- 3 x Scanning interface plate
- 1 x Blu-Tack® scan clay
- 1 x Dongle USB key
- 1 x D900 Multi-Die Fixture

### CAD Software Package

Dental System™ Premium consisting of:

- Dental Designer - Premium; see page 19  
(CAD software for dental restorations, 1 seat)
- Dental Manager  
(Order Management software, 5 seats)
- Dental Manager Inbox  
(Order reception and filtering software, 10 connections)
- TRIOS® Inbox  
(TRIOS® order reception and filtering software)

### FREE Trial Software Modules

- Abutment Designer / Implant Bars and Bridge Design  
(Limited saves: 20/5)
- Removable Partial Design  
(Limited saves: 15)
- Model Builder™ - TRIOS®  
(Limited saves: 20)

### Help System

- Integrated online user manual  
(12 languages)

## Optional:

### Add-on Modules and Module Upgrades

#### CAD Points

#### Additional Seats

#### Extended Warranty

Various add-on modules and Module Upgrades are available. See page 11-13.  
Pay-Per-Design option for selected add-on modules. Purchase in webshop.  
Additional seats can be purchased. See separate list.  
Extend the 1-year limited warranty on 3Shape hardware by 1 or 2 additional years.

### Notes:

11. Cannot be resold for a period of 12 months. Limited to 200 saves yearly.

12. Licenses are renewed upon request every year. All upgrades are included for a period of 3 years. The name of the school or university and name of the contact person are required at order placement.



# 3Shape D900L

D900L 3D scanner with Dental System™ Premium



Item Number

**20001500**

20001505 - Demo system (see note 13)

20001510 - Educational system (see note 14)

Send your order to [order@3shape.com](mailto:order@3shape.com)

## System Configuration

### 3D Scanner

D900L 3D Scanner (4 x 5 MP cameras) optimized for dental applications. With Multi-Die Scanning and Color Texture Scanning. Expanded interior space for mounted models etc.

### Desktop PC

- High performance PC (model subject to change)
- Operating system and CAD software pre-installed
- Mouse and keyboard included (monitor NOT included)

### Operating System

- Windows 8 Pro, 64-bit, English (Support for 35 languages)

### Repair Disk

- Windows 8 Pro repair disk

### Accessories

- 1 x Scanner calibration object
- 3 x Scanning interface plate
- 1 x Blu-Tack® scan clay
- 1 x Dongle USB key
- 1 x D900 Multi-Die Fixture
- 1 x Scan Height Adapter
- 1 x Occlusion Transfer Calibration Object

### CAD Software Package

Dental System™ Premium consisting of:

- Dental Designer - Premium; see page 19 (CAD software for dental restorations, 1 seat)
- Dental Manager (Order Management software, 5 seats)
- Dental Manager Inbox (Order reception and filtering software, 10 connections)
- TRIOS® Inbox (TRIOS® order reception and filtering software)

### FREE Trial Software Modules

- Abutment Designer / Implant Bars and Bridge Design (Limited saves: 20/5)
- Removable Partial Design (Limited saves: 15)
- Model Builder™ - TRIOS® (Limited saves: 20)

### Help System

- Integrated online user manual (12 languages)

### Optional:

#### Add-on Modules and Module Upgrades

#### CAD Points

#### Additional Seats

#### Extended Warranty





Various add-on modules and Module Upgrades are available. See page 11-13. Pay-Per-Design option for selected add-on modules. Purchase in webshop. Additional seats can be purchased. See separate list. Extend the 1-year limited warranty on 3Shape hardware by 1 or 2 additional years.

#### Notes:




13. Cannot be resold for a period of 12 months. Limited to 200 saves yearly.

14. Licenses are renewed upon request every year. All upgrades are included for a period of 3 years. The name of the school or university and name of the contact person are required at order placement.



## Dental System™ Subscriptions

	<b>Dental System™ Premium - Stand-alone</b> CAD software for dental restorations including Dental Manager (new dongle; see note 15).	Item number <b>80240004</b>
	<b>Dental System™ Premium - Second seat</b> CAD software for dental restorations. Additional subscription for existing dongle on a local network (same dongle; see notes 15 & 16)	Item number <b>80240003</b>
	<b>Dental System™ Standard - Stand-alone</b> CAD software for dental restorations including Dental Manager (new dongle; see note 15).	Item number <b>80240081</b>
	<b>Dental System™ Full Contour - Stand-alone</b> CAD software for dental restorations including Dental Manager (new dongle; see note 15).	Item number <b>80241200</b>

## 3Shape TRIOS® Integration

	<b>TRIOS® Inbox</b> Inbox to receive TRIOS® scans directly from dentists. Included in Dental System™ Full Contour, Standard and Premium.	Included in Dental System™
	<b>TRIOS® Inbox - Stand-alone</b> Inbox to receive TRIOS® scans. Includes STL export option and trial version of Dental System™ - TRIOS® only (20 saves).	Item number <b>80241120</b>
	<b>Dental System™ - TRIOS® only</b> Dental System™ Premium for <u>TRIOS® scans only</u> . Includes TRIOS® Inbox and Model Builder™ for TRIOS® scans.	Item number <b>80241110</b>

## 3Shape Communicate™













	<b>3Shape Communicate™</b> 3Shape Communicate™ account for online communication between lab and dentist. Includes 2 GB storage space.	Included in Dental System™
	<b>3Shape Communicate™ - additional storage space</b> 10 GB additional storage space for 3Shape Communicate™ account.	Item number <b>80249030</b>

### Notes:

15. See page 19 for a detailed description of the various software packages and system configurations.

16. Requires an active subscription of Dental Designer on the same dongle.






## Add-on CAD Modules (Require an active subscription of Dental Designer on the same dongle.)

	<b>Abutment Designer</b> CAD software for designing customized abutments and screw-retained crowns (see note 17).	Item number <b>80243010</b>
	<b>Implant Bar and Bridge Design (only for Dental System™ Premium)</b> CAD software for designing full screw-retained implant bars and bridges with direct seating on the implant (see notes 17 & 18).	Item number <b>80243020</b>
	<b>Bundle Offer</b> Abutment Designer and Implant Bar and Bridge Design bundled together in one package (see notes 17 & 19)	Item number <b>80243010 + 80243020</b>
	<b>Removable Partial Design</b> CAD software for designing removable partial dentures. Added to existing subscription (see note 17).	Item number <b>80243052</b>
	<b>Full Dentures</b> CAD software for designing Full Dentures. Added to existing subscription (see note 17).	Item number <b>80243090</b>
	<b>Implant Studio™ (Full solution) - Dental System™ Integrated - 1 separate seat</b> Add-on software package, with Implant Planner and Surgical Guide Designer. Completely integrated into Dental System™ (see notes 17, 20 & 21).	Item number <b>85240010</b>
	<b>Model Builder™ - TRIOS® - 1 separate seat</b> CAD software for designing Digital Lab Models based on TRIOS® scans and 3Shape impression scans (see notes 17 & 21).	Item number <b>80249000</b>
	<b>Model Builder™ - TRIOS® and 3<sup>rd</sup> party intra-oral scans - 1 separate seat</b> CAD software for designing digital lab models based on TRIOS® scans, Sirona Cerec scans, iTero scans, 3M True Definition and 3Shape impression scans (see notes 17 & 21).	Item number <b>80249010</b>
	<b>Model Builder™ - STL - 1 separate seat</b> CAD software for designing digital lab models based on TRIOS® scans, Sirona Cerec scans, iTero scans, 3M True Definition and 3Shape impression scans and STL files (see notes 17 & 21).	Item number <b>80249020</b>
	<b>Ortho Analyzer - 1 separate seat</b> ScanIt Orthodontics, ScanIt Ortho Impression and OrthoAnalyzer. Scanner fixtures must be purchased separately (see notes 17 & 21).	Item number <b>80242150</b>
	<b>Appliance Designer - 1 separate seat</b> CAD software for orthodontic appliances (see notes 17 & 21).	Item number <b>80245055</b>
	<b>ScanIt Impression - Unlimited seats</b> Scanning software optimized for impression scanning with Model Builder for impression scans. Impression fixture included (see notes 17 & 21).	Item number <b>80244010</b>

### Notes:

17. See page 19 for a detailed description of the various software packages and system configurations.
18. Requires Abutment Designer add-on module.
19. To obtain the Bundle Offer, both modules must be ordered simultaneously.
20. Implant Studio™ is not yet available in all countries. For more information about availability in your country, please contact 3Shape or your local 3Shape representative. Implant Studio™ is developed and manufactured by 3Shape Medical A/S.
21. Makes use of 1 Dental Designer seat during design.

## Stand-alone Modules (Not added to existing subscription - can both be added to an existing or new dongle.)







	<b>Removable Partial Design - Stand-alone</b> CAD software for designing removable partial dentures (see note 22).	Item number <b>80243051</b>
	<b>Full Dentures - Stand-alone</b> CAD software for designing Full Dentures (see note 22).	Item number <b>80243091</b>
	<b>Implant Studio™ (Full solution) - Stand-alone</b> Stand-alone software package with Implant Planner and Surgical Guide Designer. Optionally connected to Dental System™ via 3Shape Communicate™ (see notes 22 & 23).	Item number <b>85240020</b>
	<b>Implant Studio™ (Planner) - Stand-alone</b> Stand-alone software package with Implant Planning module only. Optionally connected to 'Implant Studio™ Full version' or Dental System™ via 3Shape Communicate™ (see notes 22 & 23).	Item number <b>85240030</b>
	<b>Ortho Analyzer - Stand-alone</b> ScanIt Orthodontics, ScanIt Ortho Impression and OrthoAnalyzer. Scanner fixtures must be purchased separately (see note 22).	Item number <b>80242155</b>

### Notes:




22. See page 19 for a detailed description of the various software packages and system configurations.

23. Implant Studio™ is not yet available in all countries. For more information about availability in your country, please contact 3Shape or your local 3Shape representative. Implant Studio™ is developed and manufactured by 3Shape Medical A/S.

## Module Upgrades

	<b>Sirona Scan Import/Export interface</b> Import of Sirona Cerec/Inlab scans into 3Shape Dental System™. Export of 3Shape scans to Sirona Cerec/Inlab CAD.	Item number <b>80247010 + 80247011</b>
	<b>Upgrade - Dental System™ Full Contour to Premium</b> Upgrade Dental System™ Full Contour to Dental System™ Premium. See page 19 for comparison between software packages.	Item number <b>80240088</b>
	<b>Upgrade - Dental System™ Standard to Premium</b> Upgrade Dental System™ Standard to Dental System™ Premium. See page 19 for comparison between software packages.	Item number <b>80240085</b>
 10	<b>Dental Manager +10 seats</b> Order Management software, for order handling, archiving and searching; same dongle (see note 24).	Item number <b>80241008</b>
 30	<b>Dental Manager Inbox, 30 connections</b> Order reception and filtering software, 30 connections (see note 25).	Item number <b>80241030</b>
 100	<b>Dental Manager Inbox, 100 connections</b> Order reception and filtering software, 100 connections (see note 25).	Item number <b>80241040</b>

## CAD Points













 100	<b>100 CAD Points</b> CAD Points for design of Customized Abutments, Implant Bars & Bridges, and Removable Partial Dentures.	Item number <b>80248100</b>
 200	<b>200 CAD Points</b> CAD Points for design of Customized Abutments, Implant Bars & Bridges, and Removable Partial Dentures.	Item number <b>80248101</b>
 500	<b>500 CAD Points</b> CAD Points for design of Customized Abutments, Implant Bars & Bridges, and Removable Partial Dentures.	Item number <b>80248102</b>

### Notes:

24. Requires an active subscription of Dental Designer on the same dongle.




25. When a larger Inbox is required, the full price of the new Inbox is charged and the price for the previously purchased Inbox is credited.  
 The subscription cost of the new Inbox applies from this point.

## Accessories

	<b>Scanner Interface Plates</b> Scan Interface plates for D500 and D700/800 series scanners.	Item number <b>06070009</b>
	<b>Dental 2-cast Fixture</b> Fixture for the D500 and D700/800 series scanners. Used for fixing the occlusal alignment of upper and lower cast during scanning.	Item number <b>80200002</b>
	<b>Occlusion Setup Tool incl. 4 high-precision transfer plates</b> Allows user to set up upper and lower arches prior to scanning and allowing user to skip the bite scan. Requires Dental System™ 2.9.9.2 or newer software versions. Works with all desktop scanner models.	Item number <b>80600104</b>
	<b>Dental Wax-up bridge Fixture</b> Fixture for the D500 and D700/800 series scanners. Used for fixing complete wax-up bridges in the scanner.	Item number <b>80200004</b>
	<b>Impression Fixture</b> Fixture for the D500 and D700/800 series scanners. Used for fixing dental impressions in the scanner.	Item number <b>80200005</b>
	<b>Implant Cast Fixture</b> Fixture for the D500 and D700/800 series scanners. Used for fixing implant models in the scanner.	Item number <b>80600000</b>
	<b>Blu-Tack®</b> 1 package of Blu-Tack® scan clay (90 gram).	Item number <b>02050005</b>
	<b>Pontic scanning Fixture</b> Pins for fixing of pontic/crown.	Item number <b>02070023</b>
	<b>Personal Computer</b> High performance PC (model subject to change). Mouse and keyboard included (monitor NOT included).	Item number <b>Subject to change</b>
	<b>SpaceExplorer 3D Mouse</b> SpaceExplorer 3D Mouse (3D Connexion) with 3Shape logo Six-degrees-of-freedom technology, QuickView navigation keys.	Item number <b>02019104</b>
	<b>Demonstration Suitcase</b> Suitcase for easy transportation of 3Shape scanners.	Item number <b>80000012</b>
	<b>3Shape Articulator</b> 3Shape Articulator used for articulation of digital lab models designed in 3Shape's Model Builder™.	Item number <b>80610005</b>



## Post & Core Scan Posts












Item number	Item	Compatible with	Scan post length	Scan post diameter
<b>APOL compatible</b>				
80610010	APL (6x3 scan posts included)	APOL system	Various	Various
80610011	AL1.4 (3 scan posts included)	APOL system	16 mm	1,4 mm
80610012	AS1.4 (3 scan posts included)	APOL system	12 mm	1,4 mm
80610013	AL1.6 (3 scan posts included)	APOL system	16 mm	1,6 mm
80610014	AS1.6 (3 scan posts included)	APOL system	12 mm	1,6 mm
80610015	AL1.7 (3 scan posts included)	APOL system	16 mm	1,7 mm
80610016	AS1.7 (3 scan posts included)	APOL system	12 mm	1,7 mm
<b>PIVOMATIC compatible</b>				
80610020	PIV (12x3 scan posts included)	PIVOMATIC system	Various	Various
80610021	PL1.2 (3 scan posts included)	PIVOMATIC system	16 mm	1,2 mm
80610022	PS1.2 (3 scan posts included)	PIVOMATIC system	12 mm	1,2 mm
80610023	PL1.3 (3 scan posts included)	PIVOMATIC system	16 mm	1,3 mm
80610024	PS1.3 (3 scan posts included)	PIVOMATIC system	12 mm	1,3 mm
80610025	PL1.4 (3 scan posts included)	PIVOMATIC system	16 mm	1,4 mm
80610026	PS1.4 (3 scan posts included)	PIVOMATIC system	12 mm	1,4 mm
80610027	PS1.5 (3 scan posts included)	PIVOMATIC system	16 mm	1,5 mm
80610028	PS1.5 (3 scan posts included)	PIVOMATIC system	12 mm	1,5 mm
80610029	PL1.6 (3 scan posts included)	PIVOMATIC system	16 mm	1,6 mm
80610030	PS1.6 (3 scan posts included)	PIVOMATIC system	12 mm	1,6 mm
80610031	PL1.7 (3 scan posts included)	PIVOMATIC system	16 mm	1,7 mm
80610032	PS1.7 (3 scan posts included)	PIVOMATIC system	12 mm	1,7 mm
<b>CYLINDRICAL compatible</b>				
80610040	CYL (14x3 scan posts included)	CYLINDRICAL systems*	Various	Various
80610041	CS0.90 (3 scan posts included)	CYLINDRICAL systems*	12 mm	0,9 mm
80610048	CL0.90 (3 scan posts included)	CYLINDRICAL systems*	16 mm	0,9 mm
80610042	CS1.00 (3 scan posts included)	CYLINDRICAL systems*	12 mm	1,0 mm
80610049	CL1.00 (3 scan posts included)	CYLINDRICAL systems*	16 mm	1,0 mm
80610043	CS1.14 (3 scan posts included)	CYLINDRICAL systems*	12 mm	1,14 mm
80610050	CL1.14 (3 scan posts included)	CYLINDRICAL systems*	16 mm	1,14 mm
80610044	CS1.25 (3 scan posts included)	CYLINDRICAL systems*	12 mm	1,25 mm
80610051	CL1.25 (3 scan posts included)	CYLINDRICAL systems*	16 mm	1,25 mm
80610045	CS1.40 (3 scan posts included)	CYLINDRICAL systems*	12 mm	1,4 mm
80610052	CL1.40 (3 scan posts included)	CYLINDRICAL systems*	16 mm	1,4 mm
80610046	CS1.50 (3 scan posts included)	CYLINDRICAL systems*	12 mm	1,5 mm
80610053	CL1.50 (3 scan posts included)	CYLINDRICAL systems*	16 mm	1,5 mm
80610047	CS1.75 (3 scan posts included)	CYLINDRICAL systems*	12 mm	1,75 mm
80610054	CL1.75 (3 scan posts included)	CYLINDRICAL systems*	16 mm	1,75 mm

\* E.g. ParaPost™ by Coltène Whaledent.

### Notes:

26. Post & Core Scan Posts are being sold through TRIOS A/S.

## Transfer Plates and more

	<b>Occlusion Transfer Calibration Object</b> For calibrating the transfer of position from physical articulator to virtual articulator. Will also enable you to skip the bite scan. Can be used with any transfer plate – user mounts object in gypsum to transfer position.	Item number <b>80600109</b>
	<b>Artex Transfer Plate</b> Used to transfer model position from the physical articulator to the dynamic virtual articulator in Dental System™.	Item number <b>80600100</b>
	<b>Artex Transfer Plate – D900L</b> Used to transfer model position from the physical articulator to the dynamic virtual articulator in Dental System™ (see note 27).	Item number <b>80600105</b>
	<b>Whip Mix Transfer Plate</b> Used to transfer model position from the physical articulator to the dynamic virtual articulator in Dental System™.	Item number <b>80600101</b>
	<b>Whip Mix Transfer Plate – D900L</b> Used to transfer model position from the physical articulator to the dynamic virtual articulator in Dental System™ (see note 27).	Item number <b>80600106</b>
	<b>SAM Transfer Plate</b> Used to transfer model position from the physical articulator to the dynamic virtual articulator in Dental System™.	Item number <b>80600102</b>
	<b>SAM Transfer Plate – D900L</b> Used to transfer model position from the physical articulator to the dynamic virtual articulator in Dental System™ (see note 27).	Item number <b>80600107</b>
	<b>KaVo Protar Transfer Plate</b> Used to transfer model position from the physical articulator to the dynamic virtual articulator in Dental System™.	Item number <b>80600103</b>
	<b>KaVo Protar Transfer Plate – D900L</b> Used to transfer model position from the physical articulator to the dynamic virtual articulator in Dental System™ (see note 27).	Item number <b>80600108</b>

### Notes:

27. Compatible with D900L scanner. Not compatible with scanners D500, D700, D710, D800, D810, D900.

## Accuracy Validation Objects



### ISO 12836 Objects

2 objects according to ISO 12836 standard, i.e. inlay object and bridge object. To be used as reference objects. Measured individually on CMM machine. Delivered with measuring report.

Item number  
**80600110**



### Implant Bar Validation Object

High precision implant bar object including scan flags, to control accuracy of scanner for implant bars and bridges. To be used as reference object. Measured individually on CMM machine. Delivered with measuring report.

Item number  
**80600002**

## Extended Warranty (hardware)



### Extended Warranty +1 Year

Extend the 1-year limited warranty on 3Shape hardware by 1 additional year (see note 28).

Item number  
**99000020**



### Extended Warranty +2 Years

Extend the 1-year limited warranty on 3Shape hardware by 2 more years (see note 28).

Item number  
**99000030**

### Notes:

28. Please refer to 3Shape Terms and Conditions for details.

## 3D Dental Scanners



### D500 series

Building on 3Shape's proven scanning technologies, the D500 provides a cost-efficient solution and an easy entry into CAD/CAM for labs of any size and budget.

### D700 series

The D700 is 3Shape's user-acclaimed success and a workhorse of productivity. The right choice for medium-sized labs where high throughput is the key criteria or larger labs expanding their scan-stations.

### D800 series

Texture scanning and high-res 5MP cameras make the D800 scanners the right solution for medium to large labs that demand the utmost accuracy and detail levels for even advanced indications.

### D900 series

Four high-res 5MP cameras and blue LED speckle-free technology provide extraordinary speed and highly accurate color scanning. Perfect for large high volume, productive-orientated labs.

	D500 series	D700 series	D800 series	D900 series
Scanner models	D500	D700, D710	D800, D810	D900, D900L
Cameras	2 cameras 1.3 Mega Pixel	2 cameras 1.3 Mega Pixel	2 cameras 5.0 Mega Pixel	4 cameras 5.0 Mega Pixel
Light technology	Red laser	Red laser	Red laser	Blue LED
Accuracy Crown and Bridge (ISO 12836)*	10 microns	10 microns	7 microns	7 microns
Accuracy Implant bars*	-	12 microns	8 microns	8 microns
Gypsum model scanning	✓	✓	✓	✓
Impression scanning	✓ (Add-on)	✓ (Add-on)	✓ (Add-on)	✓ (Add-on)
Multi-Die scanning	-	✓ (710 models)	✓ (810 models)	✓
Expanded interior space	-	-	-	✓ (D900L models)
Texture scanning	-	-	✓ (Standard Texture)	✓ (Color Texture)
Scan time - single die	50 Sec (55 sec)	25 sec (30 sec)	25 sec (30 sec)	15 sec (19 sec)
Scan time - 3 unit bridge	160 sec (185 sec)	80 sec (100 sec)	80 sec (100 sec)	50 Sec (70 sec)
Typical software package**	Dental System™ Standard / Full Contour	Dental System™ Premium	Dental System™ Premium	Dental System™ Premium

\* Validated and documented in accordance with ISO 12836 and with specially designed implant objects (accuracy always measured on longest span)

\*\* The availability of specific configurations depends on country and reseller

## Dental System® 2014



	Dental System™ Full Contour	Dental System™ Standard	Dental System™ Premium
CAD Indications	Copings and bridge frameworks	Upgrade to Premium	✓
	Anatomical copings and bridge framework	Upgrade to Premium	✓
	Full anatomical crowns and bridges	✓	✓
	Gingiva on all indications	Upgrade to Premium	✓
	Single wax up / Wax up bridge	✓	✓
	Inlay / Onlay / Inlay bridge / Veneers	✓	✓
	Multi-layer crowns and bridges	Upgrade to Premium	✓
	Digital temporaries	✓	✓
	Virtual diagnostic wax-ups	✓	✓
	Post & core	Upgrade to Premium	✓
	Telescopes	✓	✓
	Customized abutments	Add-on	Add-on
	Implant bars and bridges	Scanning not available for D500	Add-on
	Full dentures	Add-on	Add-on
	Removable partial dentures	Add-on	Add-on
	Customized impression trays	Upgrade to Premium	✓
	Model Builder	Add-on	Add-on
	Splints, night guards and protectors	Upgrade to Premium	✓
Advanced Tools	Orthodontics / Orthodontic appliances	Add-on	Add-on
	Implant Planning and Surgical guides (see note 29)	Add-on	Add-on
	Advanced free-form modeling	✓	✓
	Mirroring teeth	✓	✓
	Simultaneous modeling on upper & lower	✓	✓
	Dynamic virtual articulation	✓	✓
	Attachments	Upgrade to Premium	✓
	RealView™ Engine with 2D Design guidelines	✓	✓
	Align and re-use design	✓	✓
	CAD block placement	✓	✓
	3Shape Communicate™ (2 GB storage included)	✓	✓
	TRIOS® Inbox	✓	✓
	Basic / Expert mode	✓	✓
	Advanced design validation tools	✓	✓
	Manufacturing Inbox™	5 connections	10 connections
	Full manufacturing equipment compatibility	✓	✓
	3Shape LABcare™ included in subscription	✓	✓

✓	Included in the specific Dental System™ Software version
Add-on	Requires purchase of add-on module (sold separately)
Not available	Not available for this scanner model / software version
Upgrade to ...	Requires upgrade to Dental System™ Standard or Dental System™ Premium edition (sold separately)

### Notes:

29. Implant Studio™ is not yet available in all countries. For more information about availability in your country, please contact 3Shape or your local 3Shape representative. Implant Studio™ is developed and manufactured by 3Shape Medical A/S.

This Product Catalog includes products  
that are developed and manufactured by  
3Shape A/S and 3Shape Medical A/S.