

Ritter Ultimate Comfort H: left and right swivable version - the ultimate flexibility. Unit for right and left hand use.

Two versions available:

1. Left/Right, Europe version, arm mounted on water unit and
 2. Left/right, US version, arm mounted under the patient chair.
- At both versions the complete water unit can be shifted on the fly to left- and right hand treatment. Height movement range: 500-790 mm

- Right and left hand treatment versions, changeable on the fly
- Unbeatable price-ratio-performance
- Standard delivery with many features, can be equipped with camera

Standard delivery: Chairmounted, Left/Right hander version, patient chair hydraulic, synchronized, left and right armrest moveable, foot pedal, double articulated headrest, dentist element with x-ray viewer, 3 outlets for turbine or air motor without light (Midwest), 3F syringe, assistant element control panel with 3F syringe, big and small suction hoses, prepared for external suction, moveable cuspidor, water bottle system, Ritter Basic line operating lamp, 8.000 + 25.000 lux, stool MobiloFlex 40

Options: PC board for F.O outlets + one F.O. outlet (connection for max 2 F.O. outlets), 2. additional outlet with light/F.O., Borden hoses, Motor infinity brushless LED, Scaler Kolibri LED, Polylight OPM LED, Syringe Luzzani 3F, Syringe Luzzani 3F assistant, silicone pad on dentist element, up to 2 Norm trays stainless steel, sensor switch for cup filler, heating for cupfiller and instruments, Venturi Air, Place selection valve, Camera TPC (assistant side only), operating light CoolLite U or CoolLite U with sensor, stool MobiloFlex 40



Europe version



US version