



**MALMET (AUSTRALIA)
PTY LTD**

Head Office and Customer Service
ABN 95 001 717 791

9-11 McKay Avenue
PO Box 373
Leeton NSW 2705

Phone: +61 2 6953 7677
Email: info@malmet.com.au

Blanket and Fluid Warming Cabinet



Spare Parts Catalogue

Note: Due to Malmet's Policy of continuous product improvement; design and technical specifications are subject to change without notice.

Table of Contents

Layout – Front View	2
Layout – Inside Door View	3
Electrical Tray	4
Door Assembly	5
Cabinet Assembly	6
Electrical	7
Electrical Continued	8

Note: *Not all numbered parts listed in this catalogue are displayed on the Layout views.*

In the interest of continued product improvement, MALMET (AUSTRALIA) reserves the right to alter specifications without notice. E. & O. E.

Text and illustration © Copyright MALMET (AUSTRALIA) 2020

All rights reserved. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or by an information storage and retrieval system, without prior permission in writing from Malmel (Australia)


Layout – Front View



Layout – Inside Door View



Electrical Tray



Blanket & Fluid Warming Cabinets

Spare Parts Catalogue















Door Assembly

#	Part Number	Description	Diagram	#	Part Number	Description	Diagram
1	100-3015SP	Glass Door 500 X 470 PC Oyster Suitable for 105L Cabinet		1	100-3019SP	Glass Door 375 X 665 PC Oyster Suitable for 210L Cabinet	
1	83-5006SP	Glass Door 731 x 420 PC Oyster Suitable for 260L Cabinet		1	87-5001SP	Glass Door 373 X 1335 PC Oyster Suitable for 420L Cabinet	
	100-3017SP	Push-in Door Seal Gasket Suitable for 105L Cabinet			100-3018SP	Push-in Door Seal Gasket Suitable for 210L Cabinet	
	83-5008SP	Push-in Door Seal Gasket Suitable for 260L Cabinet			87-5016SP	Push-in Door Seal Gasket Suitable for 420L Cabinet	
	8-M00-83-5005SP	Door Seal Metal Strip 727mm Long			8-M00-83-5006SP	Door Seal Brush Only 727mm Long	
	100-5506SP	Door Lock Bolt Assembly					

Blanket & Fluid Warming Cabinets

Spare Parts Catalogue












Cabinet Assembly

#	Part Number	Description	Diagram	#	Part Number	Description	Diagram
	83-5035SP	PK 5 Plugs 6.35mm Plastic Smooth White			83-5040SP	PK 5 Plugs 15.88mm Plastic Smooth White	
2	83-5311SP	Inline Filter CE		3	83-5300SP	Finger Grill CE	
4	100-2235SP	Wire Shelf 414 x 455 Suitable for 105L Cabinet		5	100-2236SP	Bottom Wire Shelf 414 x 455 Suitable for 105L Fluid Warming Cabinet Only	
4	100-3000SP	Wire Shelf 655 x 455 Suitable for 210L & 420L Cabinet		5	100-3001SP	Bottom Wire Shelf 655 x 455 Suitable for 210L & 420L Fluid Warming Cabinet Only	
4	83-5075SP	Wire Shelf 414 x 455 Suitable for 260L Cabinet		5	83-5074SP	Bottom Wire Shelf 414 x 455 Suitable for 260L Fluid Warming Cabinet Only	
	83-5025SP EACH	Adjustable Feet 90mm Galv Suitable for a Freestanding Cabinet Stand			83-5030SP EACH	Cap Plastic Protective for Foot Plate Suitable for Cabinet with a Floor Skirt	
6	83-5019SP EACH	Castors – Swivel		7	83-5018SP EACH	Castors - Fixed	

Blanket & Fluid Warming Cabinets

Spare Parts Catalogue












Electrical

#	Part Number	Description	Diagram	#	Part Number	Description	Diagram
	83-5305SP	Control Temperature Thermistor CE			83-5315SP	Power Cord IEC 3 pin plug 10 Amp	
8	83-5084SP	Display Harness 9 Pin (CE – BFWF)		9	8-2050SP	Overtemp Switch (Campini 80°C) Blanket Warming Cabinet Only	
9	8-2040SP	Overtemp Switch (Campini 60°C) Fluid Warming Cabinet Only			83-5050SP	Element Heating 600 watt Round 240 volt	
10	81-2308SP	PK 5 Circuit Board Pillar (Stand Off)		11	8-K00-83-506SP	Hot Air Fan Assembly	
12	83-5002SP	Fan – 80 x 80x 25 240VAC Sunon			8-5180SP	MTS SS M3 x 10mm Phillips Head (Securing Campini Overtemp Switch)	
					92-3106SP	Screw M5x12 Pan Phillips Head 316 Ss (Securing Electrical Tray)	

Blanket & Fluid Warming Cabinets

Spare Parts Catalogue

Electrical Continued

#	Part Number	Description	Diagram	#	Part Number	Description	Diagram
13	8-M00-83-5090SP	Control Panel with BW Overlay Blanket Warming Cabinet Only		13	8-M00-84-5090SP	Control Panel with FW Overlay Fluid Warming Cabinet Only	
	83-5076SP	Display Board Blanket Warmer CE Blanket Warming Cabinet Only			84-5076SP	Display Board Fluid Warmer CE Fluid Warming Cabinet Only	
14	83-5078SP	Blanket Warmer Overlay CE Blanket Warming Cabinet Only		14	84-5078SP	Fluid Warmer Overlay CE Fluid Warming Cabinet Only	
15	83-5079SP	Relay Board (CE-Blanket/Fluid)		16	83-5083SP	Main Loom Harness CE	
	84-5095SP	Remote Sensor Harness CE		17	8-K00-84-5005SP	Remote Temp Sensor Kit Fluid Warming Cabinet Only	
					84-5096SP	Remote Temperature Sensor Fluid Warming Cabinet Only	



Manufactured by
Malmart (Australia)
ABN 95 001 717 791