


COMMITTED
TO YOUR SUCCESS.





**stryker** 

1980 81 82 83 84 85 86 87 88 89 1990 91 92 93 94 95 96 97 98 99 2000 01 02 03 04 05 06 07 08 09 2010 11 12 13 14 15

# System 7

TOTAL CONFIDENCE™

## Quality & Reliability

- Our most reliable system ever
- All handpieces must pass our 100-Point Quality Check Inspection
- No memory effect enables maximum capacity throughout the life of the battery
- System performs 100 safety checks per second when in use to ensure reliability
- Extended battery shelf life



# System 7 TOTAL CONFIDENCE™

When your team members enter the operating room, they need to feel confident that their power tools are ready to meet their needs. That's why System 7 features a newly designed Stryker SmartLiFE Battery. It runs longer, weighs less, offers built-in safety mechanisms and provides proactive warnings - so surgical teams can depend on consistency in every case.



### **ACP Technology**<sup>™</sup>



## www.www.www.Quality & Reliability

- No memory effect enables maximum capacity throughout the life of the battery
- Stryker SmartLiFE recyclable batteries provide a green solution
- Laser welded to prevent moisture from entering the battery
- Universal charger identifies overexposure to heat
- Extended battery shelf life

#### www.www.Performance

- Increased run time:
  - Small battery runs 220% longer
  - Large battery runs 183% longer
- Improved ergonomics:
  - 29% reduction in overall system weight



## **Stryker SmartLiFE™ Batteries**

- Lithium Ion cell technology
- Indicator light proactively alerts OR team when battery is near the end of its charge
- Large and small sizes compatible with all System 7 handpieces
- Stores and transfers data to universal charger to enable SORN remote device management
- Enhances performance by communicating directly with the handpiece

#### Proprietary battery safeguards provide:

- Maximum reliability by regulating voltage and alerting users in the event of overexposure in autoclave
- Unmatched safety by automatically removing voltage when the battery is not on a handpiece or charger



# System 7

TOTAL CONFIDENCE"

Real innovation is often found in small details. It's these kinds of details that make System 7 sterilization cases more cost effective, easier to use and more efficient. So you can achieve positive outcomes for your staff, your facility and your patients.

## www Efficiency

- Extended battery shelf life
- Bar codes eliminate the need for tracking labels
- Redesigned trays make it easier to load handpieces and attachments in sterilization cases
- Sealed containers minimize SPD processing time

#### Ease of Use

- Industrial grade sterilization cases are intuitively designed to be durable and easy to use
- No memory effect enables maximum capacity throughout the life of the battery

#### Cost Control

- Ability to purchase sterilization cases and inserts separately
- Stainless steel sterilization cases are more durable and last longer than previous generations

#### Education

- ProCare.com website has complete cleaning and sterilization instructions
- Cleaning and sterilization CDs available

# Automated Washer Compatible

 System designed and validated for automated washing




# System 7 Heavy Duty Power Tools

TOTAL CONFIDENCE™

# Ordering Information

System 7 brings you leading-edge technology, backed by a committed support team



# Accessories

7203-026-000	Dual Trigger Wire Collet (.7-1.8mm)
7203-126-000	Dual Trigger Pin Collet (2.0-3.2mm)
7203-226-000	Dual Trigger Pin Collet (3.0-4.2mm)
7203-036-000	Single Trigger Wire Collet (.7-1.8mm)
7203-136-000	Single Trigger Pin Collet (2.0-3.2mm)
7203-236-000	Single Trigger Pin Collet (3.0-4.2mm)
6203-110-000	A/O Quick Connect
6203-113-000	Hudson
6203-131-000	1/4" Keyed Chuck
6203-132-000	5/32" Keyed Chuck
6203-133-000	1/4" Keyless Chuck
6203-134-000	1/8" Keyless Chuck
6203-135-000	Hudson Modified Trinkle
6203-150-000	Modified Trinkle
6203-160-000	Trinkle
6203-210-000	Large A/O
6203-215-000	DHS / DCS
6203-170-000	Triathlon Adapter
7207-003-000	Sternum Blade Guard
7207-002-000	Sternum Blade Guard Long

# **Batteries**

7110-120-000	Universal Battery Charger
7212-000-000	Stryker SmartLiFE Small Battery
7215-000-000	Stryker SmartLiFE Large Battery

Steriliza	ition Cases	
7102-450-000	4 HP Sterilization Case Plus Insert Tray	
7102-454-000	3 HP Sterilization Case Plus Insert Tray	
7102-453-000	Sternum/Sabo Saw Sterilization Case Plus Insert Tray	
6102-455-000	Large Aseptic Sterilization Case	
7102-450-010	4 HP Insert Tray	
7102-454-010	3 HP Insert Tray	
7102-453-010	Sternum/Sabo Insert Tray	
7102-550-000	4 HP Sterilization Container w/Perforated Bottom	
7102-554-000	3 HP Sterilization Container w/Perforated Bottom	
7100-553-000	Sternum/Sabo Sterilization Container w/Perforated Bottom	

# *s*tryker

#### Reconstructive

Hips

Knees

Trauma & Extremities

Joint Preservation

Orthobiologics

#### **Medical & Surgical**

Power Tools & Surgical Accessories
Image Guided Navigation
Endoscopy & Arthroscopy
Integrated Communications
Beds, Stretchers & EMS
Sustainability Solutions

#### **Neurotechnology & Spine**

Craniomaxillofacial
Interventional Spine
Neurosurgical, Spine & ENT
Neurovascular
Spinal Implants



The information presented here is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with ™ designation are trademarks of Stryker. Products referenced with ® designation are registered trademarks of Stryker.

Literature Number: 9100-001-569 Rev. D PF/P.S.

Copyright © 2013 Stryker Printed in USA

Stryker Instruments 4100 East Milham Avenue Kalamazoo, MI 49001 USA t: 269 323 7700 f: 800 999 3811 toll free: 800 253 3210 www.stryker.com